# GeoMedi Co., Ltd.

Product Catalog - 2022




# GMS Intra-Oral Scanbody



- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co. guarantees durability which can be autoclavable and reusable.
- special coating allows to scan without any scan spray.
- provides with two different sizes. Short(7,5mm) / Long(11,5mm)
- consistent performance and visible in x-rays.
- Various Implant system available,
- Highly precise Geo library is provided for CAD s/w. (3shape, Exocad, Dental Wings)
- Customized Abutment Design with CNC or Premilled abutment is available.
- Zirconium prothesis design with Geo Link,





easy to verify implant system by laser marking code (Laser marking code) matching position provided identify the size of scan body by D-cut appreance

Clinical case 1



Clinical case 2




Clinical case 3






















# GMS Scanbody List


#### Osstem Implant TS (Compatible with Dio UF Regular)

| Connection Length | 11.5mm (Long) | 7.5mm (Short) |
|-------------------|---------------|---------------|
| Mini              | SUM GMS-SUML  | SUM GMS-SUMS  |
| Standard          | SURO          | SURO          |
|                   | GMS-SUROL     | GMS-SUROS     |




## Osstem Implant US

| Connection Length | 11.5mm (Long) | ı        | 7 <u>.</u> 5mm (Sh | nort)    |
|-------------------|---------------|----------|--------------------|----------|
| Mini              | EXM           | GMS-EXML | EXM J              | GMS-EXMS |
| Regular           | EXR           | GMS-EXRL | EXR                | GMS-EXRS |
|                   | 10.0mm (Long) |          |                    |          |
| Wide              | EXW           | GMS-EXWL | EXW L              | GMS-EXWS |



#### Osstem Implant SS

| Connection Length | 10mm (Long) |          | 7.5mm (Shor | t)      |
|-------------------|-------------|----------|-------------|---------|
| Regular           | ISR         | GMS-ISRL | ISR         | GMS-ISR |
| Wide              | ISW         | GMS-ISWL | ISW         | GMS-ISW |



#### Dentium Super—Line (Compatible with Megagen AnyOne & Dentis OneQ)

| Connection Length | 11.5mm (Long) | 7.5mm (Short) |     |          |
|-------------------|---------------|---------------|-----|----------|
| Universal         | DER           | GMS-DERL      | DER | GMS-DERS |



GMS Guide Driver

1.20 Hex

GMST-120H65

#### Megagen AnyRidge

| Connection Length | 11.5mm (Long) |           | 7,5mm (Short) |          |
|-------------------|---------------|-----------|---------------|----------|
| Universal         | MARR          |           | MARR          |          |
|                   |               | GMS-MARRL | GIV GIV       | IS-MARRS |

# Neo Implant IS

| Connection Length | 11.5mm (Long) |           | 7.5mm (Short) |      | GMS Guide Driver |                         |  |
|-------------------|---------------|-----------|---------------|------|------------------|-------------------------|--|
| Universal         | DERN          | GMS-DERNL | DERN GMS-DERI | NS F |                  | 1,20 Hex<br>GMST-120H65 |  |

#### Straumann Bone Level

| Connection Length | 11.5mm (Long) | 7.5mm (Short) | GMS Guide Driver  |
|-------------------|---------------|---------------|-------------------|
| SC                | BLS GMS-BLSL  |               |                   |
| NC                | BLM GMS-BLML  | BLM GMS-BLMS  | Torx<br>GMST-ST53 |
| RC                | BLR GMS-BLRL  | BLR GMS-BLRS  | Torx<br>GMST-ST65 |

GMS Guide Driver

GMST-ST65

GMST-ST65

GMST-NB53

Unigrip

GMST-NB65

Unigrip GMST—NB65

GMST-NB65

#### Straumann Tissue Level

| Connection Length | 10.0mm (Long) |          | 7.5mm (Sho | rt)      |   |  |
|-------------------|---------------|----------|------------|----------|---|--|
| RN                | ITR           | GMS-ITRL | ITR III    | GMS-ITRS |   |  |
| WN                | ITW           | GMS-ITWL | ITW        | GMS-ITWS |   |  |
|                   |               |          |            |          | - |  |

#### Nobel Biocare Active

| Connection Length | 11.5mm (Long) |          | 7.5mm (Short) |          | GMS Guide Driver |                      |  |
|-------------------|---------------|----------|---------------|----------|------------------|----------------------|--|
| NP                | ACN           | GMS-ACNL | ACN           | GMS-ACNS |                  | Unigrip<br>GMST—NB53 |  |
| RP                | ACR           | GMS-ACRL | ACR           | GMS-ACRS |                  | Unigrip<br>GMST—NB65 |  |
| WP                | ACW           | GMS-ACWL | ACW)          | GMS-ACWS |                  | Unigrip<br>GMST—NB65 |  |

## Nobel Biocare Replace

|                   | 1              |                |                  |  |  |
|-------------------|----------------|----------------|------------------|--|--|
| Connection Length | 11.5mm (Long)  | 7.5mm (Short)  | GMS Guide Driver |  |  |
| NP                | RPM GMS-RPML   | RPM GMS-RPMS   | Uniç<br>GMST-    |  |  |
| RP                | RPR GMS-RPRL   | RPR GMS-RPRS   | Uniq<br>GMST-    |  |  |
| WP                | RPW GMS-RPWL   | RPW GMS-RPWS   | unig<br>GMST-    |  |  |
| 60                | RP60 GMS-RP60L | RP60 GMS-RP60S | GMST-            |  |  |

#### Astra Tech Osseospeed TX

| Connection Length    | 11.5mm (Long) | 11.5mm (Short) | GMS Guide Driver      |  |  |
|----------------------|---------------|----------------|-----------------------|--|--|
| Aqua<br>(3.5 / 4.0)  | ASR GMS-ASRL  | ASR GMS-ASRS   | 1,27 Hex<br>GMST-127H |  |  |
| Lilac<br>(4.5 / 5.0) | ASW GMS-ASWL  | ASW GMS-ASWS   | 1,27 Hex<br>GMST-127H |  |  |
|                      |               |                |                       |  |  |

## Zimmer Tapered Screw-vent (Compatible with MIS seven 3.75/4.2, 5.0/6.0)


| _                         |               |           |               |           |                  |                         |
|---------------------------|---------------|-----------|---------------|-----------|------------------|-------------------------|
| Connection Length         | 11.5mm (Long) |           | 7.5mm (Short) |           | GMS Guide Driver |                         |
| 3.5<br>MIS seven 3.75/4.2 | TSV3          | GMS-TSV3L | TSV3          | GMS-TSV3S |                  | 1,27 Hex<br>GMST-127H53 |
| 4.5<br>MIS seven 5.0/6.0  | TSV4          | GMS-TSV4L | TSV4          | GMS-TSV4S |                  | 1,27 Hex<br>GMST-127H65 |
| 5.7                       | TSV5          | GMS-TSV5L | TSV5          | GMS-TSV5S |                  | 1.27 Hex<br>GMST—127H65 |
|                           |               |           |               |           |                  |                         |

#### Biomet 3i Certain

| Connection Length | 11.5mm (Long) |          | 7.5mm (Short) | GMS Guide Driver |                         |
|-------------------|---------------|----------|---------------|------------------|-------------------------|
| 3.4               | 3ICM          | GMS-ASRL | 3ICM GMS-ASRS |                  | 1,20 Hex<br>GMST-120H53 |
| 4.1 / 5.0 / 6.0   | 3ICR          | GMS-ASWL | 3ICR GMS-ASWS |                  | 1,20 Hex<br>GMST-120H65 |

#### MIS C1

| Connection Length | 11.5mm (Long) |          | 7,5mm (Short) |          |
|-------------------|---------------|----------|---------------|----------|
| 3,3               | C1M           | GMS-C1ML | C1M C1M       | GMS-C1MS |
| 3.75/4.2          | C1R           | GMS-C1RL | C1R           | GMS-C1RS |
| 5.0               | C1W           | GMS-C1WL | C1W           | GMS-C1WS |



# GOOMEDI

#### (주) 지 오 메 디

경기도 의왕시 이미로 40 (포일동, 인덕원IT밸리) D동 702~704, 709-1호 Tel: 02-588-2871 / Fax: 02-6280-2849 sales@geomedi.com www.geomedi.com

## GeoMedi Co., Ltd.

#709-1, 702~704, Unit D, Indeokwon IT Valley, 40, Imi-ro, Uiwang-si, Gyeonggi-do, Korea Tel: 82-2-588-2871 Fax: 82-2-6280-2849 sales@geomedi.com www.geomedi.com