





Implant Solutions
PREMIUM CAD/CAM PRODUCT



Implant Library System
(3Shape, exocad, Dentalwings)

LL/LS [Multibase]

CMFit [Premill, Blank]

3D-Analog

SC / GMS / GTS - Scanbody

Sirona Multibase

CCM Casting Abutment

Geo Screw Driver

We are exploring excellent dental implant solutions

Since 2009, we have been developing patient-customize abutment, and all of our dentists and dental technicians are also trying to make the fittest implant prosthesis for each individual person.

For Ideal Implant Restorations

We have a commitment to provide the best implant solution at a reasonable price, easy and fast way for everyone.

Co-Cr Abutment Wax up & Casting

CAD/CAM Abutment CAD design & CAM milling

For Digital Implant Solution


It is possible to make implant prosthesis by using a scanner, dental CAD, milling machine and 3D printer as developing our own digital solution for dental clinic and laboratory.

Contents

■ Company History / Certification -----	04
■ Dental CAD/CAM System Workflow -----	05
■ Design concept -----	06
GMS Scanbody /SC Scanbody -----	06
TiN screw / Geo LS Multibase -----	08
Geo LL Multibase -----	09
Geo CMFit / Reverse Jig -----	10
Geo Multi-Holder for CMFit-----	11
3D-Analog-----	12
Geo cad/cam prosthetic system-----	13
with digital library	
■ Compatible to Implant Systems-----	14
Osstem Implant TS -----	14
Osstem Implant US -----	15
Osstem Implant SS -----	16
Dentium Super-Line / Implantium -----	17
Dentium NR-Line -----	18
Megagen AnyRidge -----	19
Megagen AnyOne -----	20
Megagen EZ-plus -----	21
Megagen Ex-Feel [Internal] -----	22
Megagen MiNi -----	22
DIO SM -----	23
DIO UF / UFII -----	23
NeoBiotech IS -----	24
Dentis S-Clean / OneQ -----	25
Straumann Bone Level -----	26
Straumann Tissue Level -----	27
Straumann BLX -----	27
Nobel Biocare Active -----	28
Nobel Biocare Branemark -----	29
Nobel Biocare Replace -----	30
Dentsply Xive / Frialit2 -----	31
Dentsply Ankylos -----	32
Biomet 3i Certain -----	33
Biomet 3i Osseotite -----	34
Implant Direct 3.0 -----	35
Zimmer Tapered Screw-vent -----	36
Zimmer SPLINE -----	37
Astra Tech Osseospeed TX -----	38
Astra Tech Osseospeed EV -----	39
CAMLOG -----	40
Conelog -----	40
MIS SEVEN -----	41
MIS C1 -----	42
SIC ace / max -----	43
Thommen SPI -----	44
BEGO -----	45
MYTIS Arrow Implant -----	46
Kyocera POI -----	46
Biohorizons Laser-Lok 3.0 -----	47
Anthogyr AXIOM -----	48
WARANTEC Oneplant -----	49
Tekka In-Kone -----	50

Keystone PrimaConnex / Max-TL -----	50
Neudent CM -----	51
Neudent GM -----	51
BICON Implant -----	52
Medentika IPS -----	53
GC Aadvia Implant -----	53
BIOTECH KONTACT -----	54
Medentis medical ICX -----	55
B&B Dental -----	56
Bredent Medical Sky -----	56
NEOSS ProActive -----	57
Champions (R)evolution -----	57

■ MultiUnit(Multi-base) -----	57
■ GTS Transfer Scanbody -----	60
■ Sirona Multibase -----	62
■ CCM Casting Abutment -----	66
■ Geo Screw Driver -----	71
■ Instructions of CAD/CAM Abutment -----	74
■ Instructions of CCM Casting Abutment -----	75
■ Geo Library -----	77
■ Selection of Library -----	78
■ Geo Bridge-Single Library -----	79
■ Products List -----	81




Company History & Certification



- 2009 • 01 Released Geo CCM Abutment
- 2010 • 03 Established GeoMedi. Co., Ltd.
- 2011 • 12 Released Geo Link, CMFit, Library, Geo Driver
- 2012 • 10 The 1st Geo Symposium at Seoul National University
- 2013 • 08 Obtained KGMP
08 Obtained KFDA - implant, endosseous, laboratory supplies
10 Obtained KFDA - Geo Link Abutment, Geo CCM Abutment
12 Obtained KFDA - Geo Cmfit Abutment
- 2014 • 07 Obtained KFDA - Geo Tcp
09 The 2nd Geo Symposium at Seoul Trade Exhibition & Convention
- 2015 • 02 Obtained ISO13485, International Organization for Standardization' Medical devices Certification
02 Obtained CE, the European quality certification
04 Released Geo CAD/CAM system, Geo Mill
06 Certified as a venture business
07 Expansion & Relocation of corporation
09 Established research institute affiliated with a corporation
- 2016 • 03 Obtained KFDA - Geo Scanbody, Geo 3D Analog, Geo screw
06 Obtained KFDA - Geo Multibase Abutment
- 2017 • 07 Obtained CE, the European quality certification - Geo CCM Abutment, Geo CMFit Abutment
07 Obtained EN ISO13485 : 2012
11 Obtained certification of INNOBIZ - Ministry of SMEs and Startups
- 2018 • 05 Obtained CE, the European quality certification - Geo Multibase Abutment
11 Obtained EN ISO13485 : 2016 (31720001)
12 Achieved 1 Million dollars in export
- 2019 • 09 Obtained certification of MainBIZ - Ministry of SMEs and Startups
- 2020 • 02 Obtained KFDA - GTS Scanbody, Multi Screw Abutment
- 2021 • 04 Being selected as a Small-Ginat enterprises by Ministry of Employment and Labor
05 Obtained KFDA - OMEGA System
05 Obtained CE - GTS Scanbody & Geo Multiscrew Abutment
12 Certified as a venture business - Innovative Growth

Dental CAD/CAM System Workflow



Design concept

GMS Scanbody (Intra-oral)

- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co, guarantees durability which can be autoclavable and reusable.
- special coating allows to scan without any scan spray.
- provides with two different sizes. – Short(7.5mm) / Long(11.5mm)
- consistent performance and visible in x-rays.
- Various implant system available.
- Highly precise Geo library is provided for CAD s/w, (3shape, exocad, Dental Wings)
- Customized Abutment Design with CNC or Premilled abutment is available.
- Zirconia prothesis design with Geo Link.



GMS Short (7.5mm)

GMS Long (11.5mm)




Design concept

GMS Scanbody (Intra-oral)



GMS Scanbody / SC Scanbody

SC Scanbody(Model)




Dental Lab Scan Body.
Scan spray not necessary.
Semipermanent scanbody which assures accurate scan data.

Design concept


TiN Screw

Surface treatment of TiN-coating
- increased of abrasion resistance
- improved fixation of screw



Screw made of Titanium Grade ELI
Two Screw - TiN coating screw & Titanium screw
TiN coating screw - Final Setting

LS Multibase



Design concept


LL Multibase

LL Multibase height is editable on Geo Library.



NEW Item

Gingival height - 2.0mm
height - 6.0mm



Gingival height-1.3mm --> height-5.5mm
Gingival height-2.0mm --> height-6.0mm
Two different gingival height is available.

Zirconia - Screw retained restorations
Customized Zirconia abutment
Geo Multibase & Screw made of Titanium Grade ELI
Two Screw - TiN coating screw & Titanium screw



Geo LL Multibase(Link)

Design concept

Geo CMFit [Premill, Blank]



Diameter of $\varnothing 10$ & $\varnothing 14$ size available

Laser marked Geo & Implant System code

Provide with two different screw (TiN coating & Titanium)

Qualified material - Ti-6Al-4V ELI made by CARPENTER co.

CMFit [Premill, Blank] & Reverse Jig

Reverse Jig



Occlusal milling is possible depend on machine flexibility


Support pin is not necessary

Fix pre-milled abutment with ordinary screw

Addtional Geo holder(jig) is necessary


Design concept

Geo Multi-Holder for CMFit




Multi-holder and Reverse holder

Geo Multi-Holder for CMFit




Customized Geo Multi holder compatible with leading Milling Machines.
Available for all type of implant system.



Design concept

Geo 3D Analog

3D Analog



Pioneer of providing 3D digital Analog

Can be used both stone model & 3D printed model

Reasonable price with excellent durability

Provided all type of implant system


Stainless Steel



Cross section
3D Analog / LL -Link / Zirconia



Geo CMFit for Titan blank



Geo Link for Zirconia

Design concept

Geo cad/cam prosthetic system with digital library

Dental CAD s/w

3shape ▶

exocad

dental wings


Dental CAM s/w

FOLLOW-ME!
TECHNOLOGY GROUP
hyperDENT
COMPACT CLASSIC OPTIONS

sum3D
MILL BOX
One step beyond

vhf

Geo Scanbody




Scan body that provide the precise library matching according to the user's environment.

CMFit Abutment




Geamedi Library provides prosthesis definite sitting of a custom abutment

Multibase Abutment





Dedicated libraries for the various gingival height of Geo Ti-base abutment is well known for their emphasis on aesthetic and functional design.

3D Analog



Retention groove on 3D analog and extra space to fill up wax provides stable fix

Multifunctional analog for ordinary stone model and 3D printed model.



Precise direct interface milling from Geo geometry in cooperation with dental CAM s/w companies.

Geo ca/cam prosthetic system with digital library

Compatible to Implant Systems

Osstem Implant TS

Osstem Implant TS

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex
				20Ncm Ø3.0 / 3.5
				30Ncm Ø4.0 / 4.5 / 5.0

Description-Library code	0101 - SUM	0102 - SURO
Product Connection Platform Ø	TS Mini Ø3.0 / 3.5	TS Standard Ø4.0/4.5/5.0
SC-Scanbody	SC-SUM	SC-SURO
GMS-Scanbody(Short)	GMS-SUMS	GMS-SUROS
GMS-Scanbody(Long)	GMS-SUML	GMS-SUROL
LL-Multibase (Single)	GH : 1.3 LL2-SUM13-H2 GH : 2.0 LL2-SUM20-H2 GH : 3.0 -	LL2-SURO13-H2 LL2-SURO20-H2 LL2-SURO30-H2
LL-Multibase (Bridge)	GH : 1.3 LL2-SUM13-N2 GH : 2.0 LL2-SUM20-N2 GH : 3.0 -	LL2-SURO13-N2 LL2-SURO20-N2 LL2-SURO30-N2
LS-Multibase / Single	LS2-SUM-H2	LS2-SURO-H2
LS-Multibase / Bridge	LS2-SUM-N2	LS2-SURO-N2
CMFit Abutment(S) / Ø10	CMFit10-SUM-H2	CMFit10-SURO-H2
CMFit Abutment(B) / Ø10	CMFit10-SUM-N2	CMFit10-SURO-N2
CMFit Abutment(S) / Ø14	CMFit14-SUM-H2	CMFit14-SURO-H2
CMFit Abutment(B) / Ø14	-	CMFit14-SURO-N2
3D-Analog	3D-SUM	3D-SURO
Screw(5ea)	AC-SUM5	AC-SURO5

Compatible to Implant Systems

Osstem Implant US

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.2Hex 30Ncm

Description-Library code	0103 - EXM	0104 - EXR	0105 - EXW
Product \ Connection Platform Ø	US Mini	US Regular	US Wide
SC-Scanbody	SC-EXM	SC-EXR	SC-EXW
GMS-Scanbody(Short)	GMS-EXMS	GMS-EXRS	GMS-EXWS
GMS-Scanbody(Long)	GMS-EXML	GMS-EXRL	GMS-EXWL
LL-Multibase / Single	LL2-EXM-H2	LL2-EXR-H2	LL2-EXW-H2
LL-Multibase / Bridge	LL2-EXM-N2	LL2-EXR-N2	LL2-EXW-N2
LS-Multibase / Single	LS2-EXM-H2	LS2-EXR-H2	LS2-EXW-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXM-H2	CMFit10-EXR-H2	CMFit10-EXW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-
CMFit Abutment(S) / Ø14	CMFit14-EXM-H2	CMFit14-EXR-H2	CMFit14-EXW-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXM5	AC-EXR5	AC-EXW5

Compatible to Implant Systems

Osstem Implant SS

Osstem Implant SS

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(10mm) (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex 30Nm

Description-Library code	0106 - ISR	0107 - ISW
Connection Platform Ø	SS Regular	SS Wide
Product	Ø4.8	Ø6.0
SC-Scanbody	SC-ISRW	SC-ISRW
GMS-Scanbody(Short)	GMS-ISRS	GMS-ISWS
GMS-Scanbody(Long)	GMS-ISRL	GMS-ISWL
LL-Multibase / Single	LL2-ISR-02	LL2-ISW-02
LL-Multibase / Bridge	LL2-ISR-N2	LL2-ISW-N2
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-ISR-02	CMFit10-ISW-02
CMFit Abutment / Ø14	CMFit14-ISR-02	CMFit14-ISW-02
3D-Analog	3D-ISR	3D-ISW
Screw(5ea)	AC-EXR5	AC-EXR5

Compatible to Implant Systems

Dentium Super-Line / Implantium

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.27Hex 25~30Ncm

Description-Library code	0201 - DER	
Product	Connection Platform Ø	
SC-Scanbody	Universal	Ø3.7/4.0/4.5
GMS-Scanbody(Short)		GMS-DERS
GMS-Scanbody(Long)		GMS-DERL
LL-Multibase (Single)	GH : 1.3	LL2-DER13-H2
	GH : 2.0	LL2-DER20-H2
	GH : 3.0	LL2-DER30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-DER13-N2
	GH : 2.0	LL2-DER20-N2
	GH : 3.0	LL2-DER30-N2
LS-Multibase / Single		LS2-DER-H2
LS-Multibase / Bridge		LS2-DER-N2
CMFit Abutment(S) / Ø10		CMFit10-DER-H2
CMFit Abutment(B) / Ø10		CMFit10-DER-N2
CMFit Abutment(S) / Ø14		CMFit14-DER-H2
CMFit Abutment(B) / Ø14		CMFit14-DER-N2
3D-Analog		3D-DER
Screw(5ea)		AC-DERS

Compatible to Implant Systems

Dentium NR-Line

*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
	-	-	-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-	-		1.0square
				20Ncm

Description-Library code	0202 - DNR
Product	Connection Platform Ø Standard
OS-Scanbody	OS-DNR
LL-Multibase / Single	-
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-DNR-0
CMFit Abutment / Ø14	-
3D-Analog	-
Screw(Sea)	AC-DNR5

Compatible to Implant Systems

Megagen AnyRidge

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex 35Ncm

Description-Library code	0301 - MARR	
Product	Connection Platform Ø	Universal
		Ø3.5/4.0/4.5
SC-Scanbody	SC-MARR	
GMS-Scanbody(Short)	GMS-MARRS	
GMS-Scanbody(Long)	GMS-MARRL	
LL-Multibase (Single)	GH : 1.3	LL2-MARR-H2
	GH : 2.0	LL2-MARR20-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-MARR-N2
	GH : 2.0	LL2-MARR20-N2
LS-Multibase / Single	LS2-MARR-H2	
CMFit Abutment / Ø10	CMFit10-MARR-H2	
CMFit Abutment / Ø14	CMFit14-MARR-H2	
3D-Analog	3D-MARR	
Screw(5ea)	AC-MARR5	

Compatible to Implant Systems

Megagen AnyOne

Megagen AnyOne

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral) *Using 1.27Hex Driver*	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				25Ncm

Description-Library code	0302 - DERA	
Product	Connection Platform Ø	Universal Ø 3.5/4.0/4.5/5.0
SC-Scanbody		SC-DER
GMS-Scanbody(Short)		GMS-DERS
GMS-Scanbody(Long)		GMS-DERL
LL-Multibase (Single)	GH : 1.3	LL2-DERA13-H2
	GH : 2.0	LL2-DERA20-H2
	GH : 3.0	LL2-DERA30-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-DERA13-N2
	GH : 2.0	LL2-DERA20-N2
	GH : 3.0	LL2-DERA30-N2
LS-Multibase / Single		-
LS-Multibase / Bridge		-
CMFit Abutment(S) / Ø10		CMFit10-DERA-H2
CMFit Abutment(B) / Ø10		CMFit10-DERA-N2
CMFit Abutment(S) / Ø14		CMFit14-DERA-H2
CMFit Abutment(B) / Ø14		CMFit14-DERA-N2
3D-Analog		3D-DER
Screw(5ea)		AC-MA05

Compatible to Implant Systems

Megagen EZ-plus

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.2Hex 35Nm
	-	-		

Description-Library code	0303 - MEPM	0304 - MEPR
Connection Platform Ø	Ø 3.3	Ø 4.0/4.5/5.0
Product		
SC-Scanbody	SC-MEPM	SC-MEPR
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	-	-
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-MEPM-E	CMFit10-MEPR-E
CMFit Abutment / Ø14	-	-
3D-Analog	-	-
Screw(5ea)	AC-MEPM5	AC-MEPR5

Megagen EZ - plus

Compatible to Implant Systems

Megagen ExFeel Internal

*set (Scanbody/Multibase/CMFit Abutment + screw)

Megagen ExFeel Internal

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex
				30Nm

Description-Library code	0305 - ITRK
Connection Platform Ø	Regular
Product	Ø4.8
SC-Scanbody	SC-ITRW
GMS-Scanbody(Short)	-
GMS-Scanbody(Long)	-
LL-Multibase / Single	LL2-ITRK-02
LL-Multibase / Bridge	LL2-ITRK-N2
LS-Multibase / Single	-
CMFit Abutment(S) / Ø10	CMFit10-ITRK-02
CMFit Abutment(B) / Ø10	CMFit10-ITRK-N2
CMFit Abutment(S) / Ø14	CMFit14-ITRK-02
3D-Analog	3D-ITR
Screw(5ea)	AC-EXR5

Megagen MiNi

*set (Scanbody/Multibase/CMFit Abutment + screw)

Megagen MiNi

SC-Scanbody(Model)	Single LL-Multibase	Ø10 CMFit Abutment	Screw	Driver (Nm)
	-			1.2Hex
				15Nm Mini 3.0/3.25

Description-Library code	0306 - MMN
Connection Platform Ø	MINI
Product	Ø3.0/3.25
SC-Scanbody	SC-MMN
LL-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-MMN-H
Screw(5ea)	AC-MMNS

Compatible to Implant Systems

DIO SM

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	Single LS-Link	Screw	Driver (Ncm)
 		-		 1.7 Torx 30Ncm

Description-Library code	0401 - DSMM	0402 - DSMR
Connection Platform Ø	Mini	Regular
Scanbody	SC-DSMM	SC-DSMR
LS-Link / Single	-	-
CMFit Abutment $\varnothing 10$	CMFit10-DSMM-T	CMFit10-DSMR-T
Screw(5ea)	AC-DSMM5	AC-DSMR5

DIO UF / UFII

*set (Scanbody/Link/CMFit Abutment + screw)				
SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	 (11.5mm) (7.5mm)			-
$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex 20Ncm $\varnothing 3.0/3.3$ 30Ncm $\varnothing 3.8/4.0/4.5/5.0$

Description-Library code		0403 - DN30	0404 - SURD
Connection Platform Ø	UFII Narrow	$\varnothing 3.0/3.3$	Standard
Product			$\varnothing 3.8/4.0/4.5/5.0$
SC-Scanbody		SC-MMN	SC-SURO
GMS-Scanbody(Short)		-	GMS-SUROS
GMS-Scanbody(Long)		-	GMS-SUROL
LL-Link / Single	GH : 1.3	-	LL2-SURD-H2
	GH : 2.0	-	LL2-SURD20-H2
LL-Link / Bridge	GH : 1.3	-	LL2-SURD-N2
	GH : 2.0	-	LL2-SURD20-N2
CMFit Abutment(S) / $\varnothing 10$		CMFit10-DN30-H	CMFit10-SURD-H2
CMFit Abutment(B) / $\varnothing 10$		-	CMFit10-SURD-N2
CMFit Abutment(S) / $\varnothing 14$		-	CMFit14-SURD-H2
CMFit Abutment(B) / $\varnothing 14$		-	CMFit14-SURD-N2
3D-Analog			3D-SURO
Screw(5ea)		AC-DN305	AC-SURD5

DIO SM

DIO UF / UFII

Compatible to Implant Systems

NeoBiotech IS

NeoBiotech IS

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.2Hex
				30Nm Ø 3.5/4.0/4.5/5.0

Description-Library code		0501 - DERN	
Product		Connection Platform Ø	Standard Ø 3.5/4.0/4.5/5.0
SC-Scanbody			SC-SURO
GMS-Scanbody(Short)			GMS-DERNS
GMS-Scanbody(Long)			GMS-DERNL
LL-Multibase (Single)	GH : 1.3		LL2-DERN-H2
	GH : 2.0		LL2-DERN20-H2
LL-Multibase (Single)	GH : 1.3		LL2-DERN-N2
	GH : 2.0		LL2-DERN20-N2
LS-Multibase / Single		-	
CMFit Abutment(S) / Ø10		CMFit10-DERN-H2	
CMFit Abutment(B) / Ø10		CMFit10-DERN-N2	
CMFit Abutment(S) / Ø14		CMFit14-DERN-H2	
CMFit Abutment(B) / Ø14		CMFit14-DERN-N2	
3D-Analog		3D-DER	
Screw(5ea)		AC-DERN5	

Compatible to Implant Systems

Dentis S-Clean / One Q

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Link
	(11.5mm) (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex 30Nm

Description-Library code	0601- DERT	
Product	Connection Platform Ø	Regular Ø 3.7/4.2/4.7/5.2/6.0
SC-Scanbody		SC-DER
GMS-Scanbody(Short)		GMS-DERS
GMS-Scanbody(Long)		GMS-DERL
LL-Multibase (Single)	GH : 1.3	LL2-DERT13-H2
	GH : 2.0	-
LL-Multibase (Bridge)	GH : 1.3	LL2-DERT13-N2
	GH : 2.0	-
LS-Multibase / Single		-
CMFit Abutment(S) / Ø10		CMFit10-DERT-H2
CMFit Abutment(S) / Ø14		CMFit14-DERT-H2
3D-Analog		3D-DER
Screw(5ea)		AC-DERT5

Compatible to Implant Systems

Straumann Bone Level

Straumann Bone Level

*set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				Torx
				35Ncm

Description-Library code		0700 - BLS	0701 - BLM	0702 - BLR
Product	Connection Platform Ø	SC Ø2.9	NC Ø3.3	RC Ø4.1/4.8
SC-Scanbody		-	SC-BLM	SC-BLR
GMS-Scanbody(Short)		-	GMS-BLMS	GMS-BLRS
GMS-Scanbody(Long)		GMS-BLSL	GMS-BLML	GMS-BLRL
LL-Multibase (Single)	GH : 1.3	LL2-BLS-R	LL2-BLM-R2	LL2-BLR-R2
	GH : 2.0	-	LL2-BLM20-R2	LL2-BLR20-R2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-BLM-N2	LL2-BLR-N2
	GH : 2.0	-	LL2-BLM20-N2	LL2-BLR20-N2
LS-Multibase / Single		-	LS2-BLM-R2	LS2-BLR-R2
CMFit Abutment / Ø10		CMFit10-BLS-R	CMFit10-BLM-R2	CMFit10-BLR-R2
CMFit Abutment / Ø14		-	CMFit14-BLM-R2	CMFit14-BLR-R2
3D-Analog		3D-BLS	3D-BLM	3D-BLR
Screw(5ea)		AC-BLS5	AC-BLM5	AC-BLR5

Compatible to Implant Systems

Straumann Tissue Level

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (10mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) Torx 35Nm
			 	

Description-Library code	0703 - IT35	0704 - ITR	0705 - ITW
Connection Platform Ø	NN  Ø3.5	RN  Ø4.8	WN  Ø 6.5
Product			
SC-Scanbody	SC-IT35	SC-ITRW	SC-ITRW
GMS-Scanbody(Short)	-	GMS-ITRS	GMS-ITWS
GMS-Scanbody(Long)	-	GMS-ITRL	GMS-ITWL
LL-Multibase (Single)	GH : 1.3 GH : 2.0	- -	LL2-ITR-02 LL2-ITR20-02
LL-Multibase (Bridge)	GH : 1.3 GH : 2.0	- -	LL2-ITR-N2 LL2-ITR20-N2
LS-Multibase / Single	LS2-IT35-0	LS2-ITR-02	LS2-ITW-02
CMFit Abutment / Ø10	CMFit10-IT35-0	CMFit10-ITR-02	CMFit10-ITW-02
CMFit Abutment / Ø14	-	CMFit14-ITR-02	CMFit14-ITW-02
3D-Analog	-	3D-ITR	3D-ITW
Screw(5ea)	AC-IT355	AC-IT5	AC-IT5

Straumann BLX

*set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Single LL-Multibase	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) Torx 35Nm
 (10mm)					

Description-Library code	0710 - BLX
Connection Platform Ø	 Ø 3.5 - 6.5
Product	
GMS-Scanbody(Long)	GMS-BLXL
LL-Multibase / Single	LL2-BLX-T
CMFit Abutment / Ø14	CMFit14-BLX-T
3D-Analog	3D-BLX
Screw(5ea)	AX-BLX5

Compatible to Implant Systems

Nobel Biocare Active

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				unigrip
				15Ncm Ø3.0
				35Ncm Ø3.5/4.3/5.0/6.0

Description-Library code		0801- AC30	0802- ACN	0803- ACR	0804- ACW
Product	Connection Platform Ø	Ø3.0	Ø3.5	Ø4.3/5.0	Ø5.5
Scanbody		OS-AC30	SC-ACN	SC-ACR	SC-ACW
GMS-Scanbody(Short)		-	GMS-ACNS	GMS-ACRS	GMS-ACWS
GMS-Scanbody(Long)		-	GMS-ACNL	GMS-ACRL	GMS-ACWL
LL-Multibase (Single)	GH : 1.3	LL2-AC30-T2	LL2-ACN-H2	LL2-ACR-H2	LL2-ACW-H2
	GH : 2.0	-	LL2-ACN20-H2	LL2-ACR20-H2	-
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ACN-N2	LL2-ACR-N2	LL2-ACW-N2
	GH : 2.0	-	LL2-ACN20-N2	LL2-ACR20-N2	-
LS-Multibase / Single		LS2-AC30-T2	LS2-ACN-H2	LS2-ACR-H2	-
CMFit Abutment / Ø10		CMFit10-AC30-T2	CMFit10-ACN-H2	CMFit10-ACR-H2	CMFit10-ACW-H2
CMFit Abutment / Ø14		-	CMFit14-ACN-H2	CMFit14-ACR-H2	CMFit14-ACW-H2
3D-Analog		3D-AC30	3D-ACN	3D-ACR	3D-ACW
Screw(5ea)		AC-AC305	AC-ACN5	AC-ACR5	AC-ACW5

Compatible to Implant Systems

Nobel Biocare Branemark

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				35Ncm

Description-Library code	0805 - EXMB	0806 - EXRB	0807 - EXWB
Connection Platform Ø	NP Ø3.5	RP Ø4.1	WP Ø5.1
Product			
SC-Scanbody	SC-EXMB	SC-EXRB	SC-EXWB
GMS-Scanbody(Short)	-	-	-
GMS-Scanbody(Long)	-	-	-
LL-Multibase / Single	LL2-EXMB-H2	LL2-EXRB-H2	LL2-EXWB-H2
LL-Multibase / Bridge	LL2-EXMB-N2	LL2-EXRB-N2	LL2-EXWB-N2
LS-Multibase / Single	LS2-EXMB-H2	LS2-EXRB-H2	LS2-EXWB-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXMB-H2	CMFit10-EXRB-H2	CMFit10-EXWB-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXRB-N2	-
CMFit Abutment(S) / Ø14	CMFit14-EXMB-H2	CMFit14-EXRB-H2	CMFit14-EXWB-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXMB5	AC-EXRB5	AC-EXWB5

Compatible to Implant Systems

Nobel Biocare Replace

Nobel Biocare Replace

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) unigrip 35Nm

Description-Library code	0808- RPM	0809- RPR	0810- RPW	0811- RP60
Product \ Connection Platform Ø	NP Ø3.5	RP Ø4.3	WP Ø5.0	Ø 6.0
Scanbody	SC-RPM	SC-RPR	SC-RPW	SC-RP60
GMS-Scanbody(Short)	GMS-RPMS	GMS-RPRS	GMS-RPWS	GMS-RP60S
GMS-Scanbody(Long)	GMS-RPML	GMS-RPRL	GMS-RPWL	GMS-RP60L
LL-Multibase / Single	LL2-RPM-E2	LL2-RPR-E2	LL2-RPW-E2	LL2-RP60-E2
LL-Multibase / Bridge	LL2-RPM-N2	LL2-RPR-N2	LL2-RPW-N2	LL2-RP60-N2
LS-Multibase / Single	LS2-RPM-E2	LS2-RPR-E2	LS2-RPW-E2	LS2-RP60-E2
CMFit Abutment / Ø10	CMFit10-RPM-E2	CMFit10-RPR-E2	CMFit10-RPW-E2	CMFit10-RP60-E2
CMFit Abutment / Ø14	CMFit14-RPM-E2	CMFit14-RPR-E2	CMFit14-RPW-E2	CMFit14-RP60-E2
3D-Analog	3D-RPM	3D-RPR	3D-RPW	3D-RP60
Screw(5ea)	AC-RPM5	AC-RPR5	AC-RPW5	AC-RP605

Compatible to Implant Systems

Dentsply Xive/Frialit2

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.22Hex 24Ncm


Description-Library code	0901 - XV3	0902 - XVM	0903 - XVR	0904 - XVW	0905 - XV55
Product \ Connection Platform Ø	Ø3.0	Ø3.4	Ø3.8	Ø4.5	Ø5.5
SC-Scanbody	SC-XV3	SC-XVM	SC-XVR	SC-XWV	SC-XWV
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-XV3-H2	LL2-XVM-H2	LL2-XVR-H2	LL2-XWV-H2	LL2-XW55-H2
LL-Multibase / Bridge	LL2-XV3-N2	LL2-XVM-N2	LL2-XVR-N2	LL2-XWV-N2	LL2-XW55-N2
LS-Multibase / Single	LS2-XV3-H2	LS2-XVM-H2	LS2-XVR-H2	LS2-XWV-H2	LS2-XW55-H2
CMFit Abutment / Ø10	CMFit10-XV3-H2	CMFit10-XVM-H2	CMFit10-XVR-H2	CMFit10-XWV-H2	CMFit10-XW55-H2
CMFit Abutment / Ø14	-	CMFit14-XVM-H2	CMFit14-XVR-H2	CMFit14-XWV-H2	CMFit14-XW55-H2
3D-Analog	3D-XV3	3D-XVM	3D-XVR	3D-XWV	3D-XW55
Screw(5ea)	AC-XV35	AC-XV5	AC-XV5	AC-XV5	AC-XV5

Compatible to Implant Systems

Dentsply Ankylos

Dentsply Ankylos

*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
			-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.0Hex 15Nm

Description-Library code	0906 - AKLL	
Connection Platform Ø	universal  Ø 3.4/4.5/5.5	
Product	OS-Scanbody	
LL-Multibase (Single)	GH : 1.3	LL2-AKLL-T
	GH : 2.0	LL2-AKLL20-T
LL-Multibase (Bridge)	GH : 1.3	LL2-AKLL-N
	GH : 2.0	LL2-AKLL20-N
LS-Multibase / Single	-	
CMFit Abutment(S) / Ø10	CMFit10-AKLL-T	
CMFit Abutment(B) / Ø10	CMFit10-AKLL-N	
3D-Analog	3D-AKLL	
Screw	-	

Compatible to Implant Systems

Biomet 3i Certain

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.2Hex 20Ncm

Description-Library code	1001 - 3ICM	1002 - 3ICR	1003 - 3ICW	1004 - 3IC60
Connection Platform Ø	Ø3.4	Ø4.1	Ø5.0	Ø6.0
Product				
SC-Scanbody	SC-3ICM	SC-3ICR	SC-3ICR	SC-3ICR
GMS-Scanbody(Short)	GMS-3ICMS	GMS-3ICRS	GMS-3ICRS	GMS-3ICRS
GMS-Scanbody(Long)	GMS-3ICML	GMS-3ICRL	GMS-3ICRL	GMS-3ICRL
LL-Multibase / Single	LL2-3ICM-H2	LL2-3ICR-H2	LL2-3ICW-H2	LL2-3IC60-H2
LL-Multibase / Bridge	LL2-3ICM-N2	LL2-3ICR-N2	LL2-3ICW-N2	LL2-3IC60-N2
LS-Multibase / Single	LS2-3ICM-H2	LS2-3ICR-H2	LS2-3ICW-H2	LS2-3IC60-H2
CMFit Abutment / Ø10	CMFit10-3ICM-H2	CMFit10-3ICR-H2	CMFit10-3ICW-H2	CMFit10-3IC60-H2
CMFit Abutment / Ø14	-	CMFit14-3ICR-H2	CMFit14-3ICW-H2	CMFit14-3IC60-H2
3D-Analog	3D-3ICM	3D-3ICR	3D-3ICW	3D-3IC60
Screw(Sea)	AC-3IC5	AC-3IC5	AC-3IC5	AC-3IC5

Compatible to Implant Systems

Biomet 3i Osseotite

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.2Hex 20Ncm

Description-Library code	1005- EX3M	1006- EXR	1007- EXT	1008- EXTW
Connection Platform Ø	Ø3.4	Ø4.1	Ø5.0	Ø6.0
Product				
Scanbody	SC-EX3M	SC-EXR	SC-EXT	SC-EXTW
GMS-Scanbody(Short)	-	GMS-EXRS	GMS-EXRS	GMS-EXRS
GMS-Scanbody(Long)	-	GMS-EXRL	GMS-EXRL	GMS-EXRL
LL-Multibase / Single	LL2-EX3M-H2	LL2-EXR-H2	LL2-EXT-H2	LL2-EXTW-H2
LL-Multibase / Bridge	LL2-EX3M-N2	LL2-EXR-N2	LL2-EXT-N2	LL2-EXTW-N2
LS-Multibase / Single	LS2-EX3M-H2	LS2-EXR-H2	LS2-EXT-H2	LS2-EXTW-H2
LS-Multibase / Bridge	-	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EX3M-H2	CMFit10-EXR-H2	CMFit10-EXT-H2	CMFit10-EXTW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-	-
CMFit Abutment(S) / Ø14	-	CMFit14-EXR-H2	CMFit14-EXT-H2	CMFit14-EXTW-H2
3D-Analog	-	3D-EXR	-	-
Screw(5ea)	AC-EXR5	AC-EXR5	AC-EXR5	AC-EXR5

Compatible to Implant Systems

Implant Direct 3.0

*set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase	
		-	-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.25Hex 30Nm

Description-Library code	1101 - TSVM
Product	Connection Platform Ø
OS-Scanbody	Implant Direct  Ø3.0
LL-Multibase / Single	OS-TSVM
LL-Multibase / Bridge	LL2-TSVM-H2
LS-Multibase / Single	-
CMFit Abutment / Ø10	-
CMFit Abutment / Ø14	CMFit10-TSVM-H2
3D-Analog	-
Screw(5ea)	3D-TSVM
	AC-TSVM5

Compatible to Implant Systems

Zimmer Tapered Screw-vent

Zimmer Tapered Screw-vent

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.25Hex 30Ncm

Description-Library code		1201- TSV3	1202- TSV4	1203- TSV5
Product	Connection Platform Ø			
Scanbody		SC-TSV3	SC-TSV4	SC-TSV5
GMS-Scanbody(Short)		GMS-TSV3S	GMS-TSV4S	GMS-TSV5S
GMS-Scanbody(Long)		GMS-TSV3L	GMS-TSV4L	GMS-TSV5L
LL-Link / Single	GH : 1.0	LL2-TSV3-H2	LL2-TSV4-H2	LL2-TSV5-H2
	GH : 2.0	LL2-TSV320-H2	LL2-TSV420-H2	-
LL-Link / Bridge	GH : 1.0	LL2-TSV3-N2	LL2-TSV4-N2	LL2-TSV5-N2
	GH : 2.0	LL2-TSV320-N2	LL2-TSV420-N2	-
LS-Link / Single		LS2-TSV3-H2	LS2-TSV4-H2	LS2-TSV5-H2
CMFit Abutment / Ø10		CMFit10-TSV3-H2	CMFit10-TSV4-H2	CMFit10-TSV5-H2
CMFit Abutment / Ø14		CMFit14-TSV3-H2	CMFit14-TSV4-H2	CMFit14-TSV5-H2
3D-Analog		3D-TSV3	3D-TSV4	3D-TSV5
Screw(5ea)		AC-TSV5	AC-TSV5	AC-TSV5

Compatible to Implant Systems

Zimmer Spline

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
-		-		 1.25Hex 30Nm

Description-Library code	1204 - SP40	1205- SP50
Product \ Connection Platform Ø	Ø 3.7 / 4.0	Ø 5.0
SC-Scanbody	SC-SP40	SC-SP40
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SP40-H	LL2-SP50-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	-	-
CMFit Abutment / Ø14	CMFit14-SP40-H	CMFit14-SP50-H
3D-Analog	-	-
Screw(5ea)	AC-TSV5	AC-TSV5

Zimmer Spline

Compatible to Implant Systems

Astra Tech Osseospeed TX

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.27Hex 15Ncm Ø3.0 20Ncm Ø3.5/4.0 25Ncm Ø4.5/5.0

Description-Library code		1301- AS3	1302- ASR	1303- ASW
Product	Connection Platform Ø	Yellow Ø3.0	Aqua Ø3.5/4.0	Lilac Ø4.5/5.0
Scanbody		SC-AS3	SC-ASR	SC-ASW
GMS-Scanbody(Short)		-	GMS-ASRS	GMS-ASWS
GMS-Scanbody(Long)		-	GMS-ASRL	GMS-ASWL
LL-Multibase (Single)	GH : 1.3	LL2-AS3-T2	LL2-ASR-H2	LL2-ASW-H2
	GH : 2.0	-	LL2-ASR20-H2	LL2-ASW20-H2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ASR-N2	LL2-ASW-N2
	GH : 2.0	-	LL2-ASR20-N2	LL2-ASW20-N2
LS-Multibase / Single		-	LS2-ASR-H2	LS2-ASW-H2
CMFit Abutment / Ø10		CMFit10-AS3-T2	CMFit10-ASR-H2	CMFit10-ASW-H2
CMFit Abutment / Ø14		-	CMFit14-ASR-H2	CMFit14-ASW-H2
3D-Analog		-	3D-ASR	3D-ASW
Screw(5ea)		AC-AS35	AC-ASR5	AC-ASW5

Compatible to Implant Systems

Astra Tech Osseospeed EV

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.27Hex 25Ncm

Description-Library code	1304 - AEV30	1305 - AEV36	1306 - AEV42	1307 - AEV48	1308 - AEV54
Connection Platform Ø	Green	Purple	Yellow	Blue	Brown
SC-Scanbody	SC-AEV30	SC-AEV36	SC-AEV42	SC-AEV48	SC-AEV54
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-AEV30-T2	LL2-AEV36-T2	LL2-AEV42-T2	LL2-AEV48-T2	LL2-AEV54-T2
LL-Multibase / Bridge	LL2-AEV30-N2	LL2-AEV36-N2	LL2-AEV42-N2	LL2-AEV48-N2	LL2-AEV54-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-AEV30-T2	CMFit10-AEV36-T2	CMFit10-AEV42-T2	CMFit10-AEV48-T2	CMFit10-AEV54-T2
CMFit Abutment / Ø14	-	CMFit14-AEV36-T2	CMFit14-AEV42-T2	CMFit14-AEV48-T2	CMFit14-AEV54-T2
3D-Analog	3D-AEV30	3D-AEV36	3D-AEV42	3D-AEV48	3D-AEV54
Screw(5ea)	AC-AEV305	AC-AEV365	AC-AEV425	AC-AEV48545	AC-AEV48545

Compatible to Implant Systems

Camlog

Camlog

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.28Hex 20Nm

Description-Library code	1401- CAM33	1402- CAM38	1403- CAM43	1404- CAM50	1405- CAM60
Product	Connection Platform Ø	Ø 3.3	Ø 3.8	Ø 4.3	Ø 5.0
Scanbody	SC-CAM33	SC-CAM38	SC-CAM43	SC-CAM5060	SC-CAM5060
LL-Multibase / Single	LL2-CAM33-T2	LL2-CAM38-T2	LL2-CAM43-T2	LL2-CAM50-T2	LL2-CAM60-T2
LL-Multibase / Bridge	LL2-CAM33-N2	LL2-CAM38-N2	LL2-CAM43-N2	LL2-CAM50-N2	LL2-CAM60-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-CAM33-T2	CMFit10-CAM38-T2	CMFit10-CAM43-T2	CMFit10-CAM50-T2	CMFit10-CAM60-T2
CMFit Abutment / Ø14	-	CMFit14-CAM38-T2	CMFit14-CAM43-T2	CMFit14-CAM50-T2	CMFit14-CAM60-T2
3D-Analog	3D-CAM33	3D-CAM38	3D-CAM43	3D-CAM50	-
Screw(5ea)	AC-CAM5	AC-CAM5	AC-CAM5	AC-CAM50605	AC-CAM50605

Conelog

Conelog

*set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
					1.28Hex 20Nm

Description-Library code		1406- CN33	1407- CN3843	1408- CN50
Product	Connection Platform Ø	Ø 3.3	Ø 3.8/4.3	Ø 5.0/6.0
GMS-Scanbody(Long)	GMS-CN33L	GMS-CN3843L	GMS-CN50L	-
CMFit Abutment / Ø10	CMFit10-CN33-E	-	-	-
CMFit Abutment / Ø14	-	CMFit14-CN3843-E	CMFit14-CN50-E	-
3D-Analog	3D-CN33	3D-CN3843	3D-CN50	-
Screw(5ea)	AC-CN5	AC-CN5	AC-CN5	AC-CN505

Compatible to Implant Systems

MIS SEVEN

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	(11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex
				30Nm

Description-Library code	1501 - MS33	1502 - MSR	1503 - MSW
Product Connection Platform Ø	Narrow Ø3.3	Standard Ø 3.75 / 4.2	Wide Ø 5.0 / 6.0
SC-Scanbody	SC-MS33	SC-TSV3	SC-TSV4
GMS-Scanbody(Short)	-	GMS-TSV3S	GMS-TSV4S
GMS-Scanbody(Long)	-	GMS-TSV3L	GMS-TSV4L
LL-Multibase / Single	LL2-MS33-H	LL2-MSR-H	LL2-MSW-H
LL-Multibase / Bridge	LL2-MS33-N	LL2-MSR-N	LL2-MSW-N
LS-Multibase / Single	LS2-MS33-H	LS2-MSR-H	LS2-MSW-H
CMFit Abutment / Ø10	CMFit10-MS33-H	CMFit10-MSR-H	CMFit10-MSW-H
CMFit Abutment / Ø14	-	CMFit14-MSR-H	CMFit14-MSW-H
3D-Analog	3D-MS33	3D-TSV3	3D-TSV4
Screw(5ea)	AC-MS33	AC-MS5	AC-MS5

MIS SEVEN

Compatible to Implant Systems

MIS C1

MIS C1

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.25Hex 30Nm

Description-Library code	1504- C1M	1505- C1R	1506- C1W
Product \ Connection Platform Ø	C1 Ø 3.3	C1 Ø 3.75/4.2	C1 Ø 5.0
Scanbody	SC-C1M	SC-C1R	SC-C1W
GMS-Scanbody(Short)	GMS-C1MS	GMS-C1RS	GMS-C1WS
GMS-Scanbody(Long)	GMS-C1ML	GMS-C1RL	GMS-C1WL
LL-Multibase / Single	LL2-C1M-T	LL2-C1R-T	LL2-C1W-T
LL-Multibase / Bridge	LL2-C1M-N	LL2-C1R-N	LL2-C1W-N
LS-Multibase / Single	-	-	-
CMFit Abutment / Ø10	CMFit10-C1M-T	CMFit10-C1R-T	CMFit10-C1W-T
CMFit Abutment / Ø14	CMFit14-C1M-T	CMFit14-C1R-T	CMFit14-C1W-T
3D-Analog	3D-C1M	3D-C1R	3D-C1W
Screw(5ea)	AC-C1M5	AC-C15	AC-C15

Compatible to Implant Systems

SIC ace / max

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Link	Bridge LL-Link	Single LS-Link
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.2Hex 20Ncm

Description-Library code	1601- SICM	1602- SICR
Connection Platform Ø	Ø 3.3	Ø 4.2
Product		
Scanbody	SC-SIC	SC-SIC
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Link / Single	LL2-SICM-H	LL2-SICR-H
LL-Link / Bridge	LL2-SICM-N	LL2-SICR-N
LS-Link / Single	-	-
CMFit Abutment / Ø10	CMFit10-SICM-H	CMFit10-SICR-H
CMFit Abutment / Ø14	CMFit14-SICM-H	CMFit14-SICR-H
3D-Analog	3D-SICM	3D-SICR
Screw(5ea)	AC-SIC5	AC-SIC5

SIC ace / max

Compatible to Implant Systems

Thommen SPI

Thommen SPI

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				Four lobe
				15Ncm Ø3.5 25Ncm Ø 4.0/4.5/5.0/6.0

Description-Library code	1701- SPI35	1702- SPI40	1703- SPI45	1704- SPI50	1705- SPI60
Product \ Connection Platform Ø	Ø3.5	Ø4.0	Ø4.5	Ø5.0	Ø6.0
Scanbody	SC-SPI35	SC-SPI40	SC-SPI45	SC-SPI45	SC-SPI45
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-SPI35-H	LL2-SPI40-H	LL2-SPI45-H	LL2-SPI50-H	LL2-SPI60-H
LL-Multibase / Bridge	LL2-SPI35-N	LL2-SPI40-N	LL2-SPI45-N	LL2-SPI50-N	LL2-SPI60-N
LS-Multibase / Single	LS2-SPI35-H	LS2-SPI40-H	LS2-SPI45-H	LS2-SPI50-H	LS2-SPI60-H
LS-Multibase / Bridge	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-SPI35-H	CMFit10-SPI40-H	CMFit10-SPI45-H	CMFit10-SPI50-H	CMFit10-SPI60-H
CMFit Abutment / Ø14	-	CMFit14-SPI40-H	CMFit14-SPI45-H	CMFit14-SPI50-H	CMFit14-SPI60-H
3D-Analog	3D-SPI35	3D-SPI40	3D-SPI45	3D-SPI50	-
Screw(5ea)	AC-SPI35	AC-SPI40	AC-SPI45	AC-SPI50	AC-SPI55

Compatible to Implant Systems

BEGO

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.25Hex 30Ncm

Description-Library code	1801 - BEGO37	1802 - BEGO41	1803 - BEGO45	1804 - BEGO55
Connection Platform Ø				
Product				
SC-Scanbody	SC-BEGO	SC-BEGO	SC-BEGO	SC-BEGO
GMS-Scanbody(Short)	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-
LL-Multibase / Single	LL2-BEG037-H	LL2-BEG041-H	LL2-BEG045-H	LL2-BEG055-H
LL-Multibase / Bridge	LL2-BEG037-N	LL2-BEG041-N	-	-
LS-Multibase / Single	LS2-BEG037-H	LS2-BEG041-H	-	-
CMFit Abutment / Ø10	CMFit10-BEG037-H	CMFit10-BEG041-H	CMFit10-BEG045-H	CMFit10-BEG055-H
CMFit Abutment / Ø14	-	-	CMFit14-BEG045-H	CMFit14-BEG055-H
3D-Analog	3D-BEG037	-	-	-
Screw(5ea)	AC-BEG05	AC-BEG05	AC-BEG05	AC-BEG05

Compatible to Implant Systems

MYTIS Arrow Implant

MYTIS Arrow Implant

*set (Scanbody/Link/CMFit Abutment + screw)

OS-Scanbody(Intra Oral)	Ø10 CMFit Abutment	Screw	Driver (Nm)
			 1.3Hex
			30Nm

Description-Library code	1901 - AR3740	1902 - AR4760
Product	Connection Platform Ø	Connection Platform Ø
OS-Scanbody	 Ø 3.7/4.0	 Ø 4.7/6.0
	OS-AR3740	OS-AR4760
	CMFit10-AR3740-H	CMFit10-AR4760-H
Screw(5ea)	AC-AR5	AC-AR5

Kyocera POI

Kyocera POI

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Link	Bridge LL-Link	Ø10 CMFit Abutment	Screw	Driver (Nm)
		-			 1.1Hex
					20Nm

Description-Library code	2001 - POI37	2002 - POI42
Product	Connection Platform Ø	Connection Platform Ø
SC-Scanbody	 Ø 3.7	 Ø 4.2
LL-Link / Single	SC-POI37	SC-POI42
LL-Link / Bridge	LL2-POI37-H	LL2-POI42-H
CMFit Abutment / Ø10	CMFit10-POI37-H	CMFit10-POI42-H
Screw(5ea)	AC-POI375	AC-POI425

Compatible to Implant Systems

Biohorizons Laser Lok 3.0

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.25Hex 30Ncm

Description-Library code	2101 - BH30
Product	Connection Platform Ø
OS-Scanbody	 Ø 3.0
LL-Multibase / Single	LL2-BH30-H
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-BH30-H
CMFit Abutment / Ø14	-
3D-Analog	3D-BH30
Screw(5ea)	AC-BH305

Compatible to Implant Systems

Anthogyr AXIOM

Anthogyr AXIOM

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.18Hex 25Nm

Description-Library code	2201 - ANG
Product	Connection Platform Ø universal  Ø 3.4/4.0/4.6/5.2
SC-Scanbody	SC-ANG
GMS-Scanbody(Short)	-
GMS-Scanbody(Long)	-
LL-Multibase / Single	LL2-ANG-E
LL-Multibase / Bridge	LL2-ANG-N
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-ANG-E
CMFit Abutment / Ø14	-
3D-Analog	3D-ANG
Screw(5ea)	AC-ANG5

Compatible to Implant Systems

Warantec Oneplant

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	Single LS-Multibase	Screw	Driver (Nm)
		-		 1.2Hex 30Nm

Description-Library code	2301 - WOP		
Product	Connection Platform Ø	Standard	 Ø 4.3/4.8/5.3
Scanbody	SC-WOP		
LS-Multibase / Single	-		
CMFit Abutment $\varnothing 10$	CMFit10-WOP-0		
Screw(5ea)	AC-WOPS		

Tekka In-Kone

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
		-	-	
		-	-	
$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	3D-Analog	Screw	 1.2Hex 15Nm $\varnothing 3.5/4.0/4.5/5.5$
				

Description-Library code	2401 - TEK		
Product	Connection Platform Ø	Standard	 Ø 3.5/4.0/4.5/5.5
SC-Scanbody	SC-TEK		
GMS-Scanbody(Short)	-		
GMS-Scanbody(Long)	-		
LS-Multibase / Single	LS2-TEK-H		
CMFit Abutment / $\varnothing 10$	CMFit10-TEK-H		
CMFit Abutment / $\varnothing 14$	-		
3D-Analog	3D-TEK		
Screw(5ea)	AC-TEK5		

Warantec Oneplant

Tekka In-Kone

Compatible to Implant Systems

Keystone PrimaConnex / Max-TL

Keystone PrimaConnex / Max-TL

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.25 Square 30Ncm
				

Description-Library code	2501 - KEP35	2502 - KEP41	2503 - KEP50	2504 - KEX50	2505 - KEX60
Connection Platform Ø	Ø 3.5	Ø 4.1	Ø 5.0 / 6.0	Ø 7.0	Ø 8.0
Product	SC-KEP35	SC-KEP41	SC-KEP50	SC-KEX50	SC-KEX60
SC-Scanbody					
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-KEP35-T	LL2-KEP41-T	LL2-KEP50-T	LL2-KEX50-T	LL2-KEX65-T
LL-Multibase / Bridge	LL2-KEP35-N	LL2-KEP41-N	LL2-KEP50-N	LL2-KEX50-N	LL2-KEX65-N
LS-Multibase / Single	LS2-KEP35-T	LS2-KEP41-T	LS2-KEP50-T	-	-
CMFit Abutment / Ø10	CMFit10-KEP35-T	CMFit10-KEP41-T	CMFit10-KEP50-T	-	-
CMFit Abutment / Ø14	CMFit14-KEP35-T	CMFit14-KEP41-T	CMFit14-KEP50-T	CMFit14-KEX50-T	CMFit14-KEX65-T
3D-Analog	3D-KEP35	3D-KEP41	3D-KEP50	3D-KEX50	3D-KEX65
Screw(5ea)	AC-KEP5	AC-KEP5	AC-KEP5	AC-KEP5	AC-KEP5

Compatible to Implant Systems

Neodent CM

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				0.9Hex 15Ncm

Description-Library code	2601 - NCM		
Product	Connection Platform Ø	Standard	Ø 3.5/4.3/5.0
SC-Scanbody	SC-NCM		
LS-Multibase / Single	-		
CMFit Abutment / Ø10	CMFit10-NCM-H		
CMFit Abutment / Ø14	-		
3D-Analog	3D-NCM		
Screw	-		

Neodent GM

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				unigrip (not compatible with Nobel) 20Ncm

Description-Library code	2602 - NCG		
Product	Connection Platform Ø	Standard	Ø 3.5/3.75/4.0/4.3/5.0
SC-Scanbody	SC-NCG		
CMFit Abutment / Ø10	CMFit10-NCG-H		
CMFit Abutment / Ø14	CMFit14-NCG-H		
LL-Multibase / Single	LL2-NCG-H		
LL-Multibase / Bridge	-		
3D-Analog	3D-NCG		
Screw(5ea)	AC-NCG5		

Neodent CM

Neodent GM

Compatible to Implant Systems

Bicon Implant

Bicon Implant


*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Single & Bridge LL-Multibase	Single & Bridge LS-Multibase	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment
		-		-
Description-Library code				
SC-Scanbody	2701- BC20	2702- BC25	2703- BC30	
LL-Multibase / Single & Bridge	<input type="radio"/> $\varnothing 2.0$	<input type="radio"/> $\varnothing 2.5$	<input type="radio"/> $\varnothing 3.0$	
LS-Multibase / Single & Bridge	-	-	-	
CMFit Abutment / $\varnothing 10$	CMFit10-BC20-N	CMFit10-BC25-N	CMFit10-BC30-N	
CMFit Abutment / $\varnothing 14$	-	-	-	

Compatible to Implant Systems

Medentika IPS

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	Screw	Driver (Nm)
		-		 1.3 Hex
				25Nm

Description-Library code	2801- MTR
Product	Connection Platform Ø  Ø3.5/4.3/5.0
SC-Scanbody	SC-MTR
CMFit Abutment / Ø10	CMFit10-MTR-R
Screw(5ea)	AC-MTR5

Medentika IPS

GC Aadvia Implant

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	$\varnothing 10$ CMFit Abutment	$\varnothing 14$ CMFit Abutment	Screw	Driver (Nm)
		-		 Torx
				20Nm

Description-Library code	2901- GAD
Product	Connection Platform Ø  Regular / Wide Ø4.0 / 5.0
SC-Scanbody	SC-GAD
CMFit Abutment / Ø10	CMFit10-GAD-H
Screw(5ea)	AC-GAD5

GC Aadvia Implant

Compatible to Implant Systems

BIOTECH KONTACT

*set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
	-			 1.2 Hex 20Ncm

Description-Library code	3001- BK30	3002- BK40
Product	Connection Platform Ø	
SC-Scanbody	 Ø 3.0	 Ø 3.6 / 4.2 / 4.8 / 5.4
LL-Multibase / Single	LL2-BK30-T	LL2-BK40-T
LL-Multibase / Bridge	-	-
LS-Multibase / Single	LS2-BK30-T	LS2-BK40-T
CMFit Abutment / Ø10	CMFit10-BK30-T	CMFit10-BK40-T
3D-Analog	3D-BK30	3D-BK40
Screw(5ea)	AC-BK305	AC-BK405

Compatible to Implant Systems

Medentis medical ICX-System

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm)
				1.42Hex 30Nm

Description-Library code	3101 - ICX
Product	Connection Platform Ø
SC-Scanbody	universal  Ø3.75/4.1/4.8
LL-Multibase / Single	LL2-ICX-H
LL-Multibase / Bridge	LL2-ICX-N
LS-Multibase / Single	LS2-ICX-H
CMFit Abutment / Ø10	CMFit10-ICX-H
CMFit Abutment / Ø14	CMFit14-ICX-H
3D-Analog	3D-ICX
Screw(5ea)	AC-ICX5

Compatible to Implant Systems


B&B Dental

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-			-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.27 Hex 25Nm

Description-Library code	3201- BNB
Connection Platform Ø	 Ø 4.0/4.5/5.0
Product	
GMS-Scanbody(long)	GMS-BNBL
LL-Multibase / Single	LL2-BNB-H
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	CMFit10-BNB-H
3D-Analog	3D-BNB
Screw(5ea)	AC-BNB5

Bredent Medical Sky

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-			-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx 25Nm

Description-Library code	3301- SKR
Connection Platform Ø	 Ø 3.5(NP) /4.0-6.0
Product	
GMS-Scanbody(long)	GMS-SKRL
LL-Multibase / Single	LL2-SKR-T
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	CMFit10-SKR-T
CMFit Abutment / Ø14	CMFit14-SKR-T
3D-Analog	3D-SKR
Screw(5ea)	AC-SKR5

B&B Dental

Bredent Medical Sky

56

Compatible to Implant Systems

NEOSS ProActive

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase			
-		-	-	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) (not compatible with Nobel)			
-				32Nm			
Description-Library code				3401- NEP			
Product	Connection Platform Ø	Ø 3.5/4.0/4.5/5.0/5.5					
GMS-Scanbody(long)	GMS-NEPL						
LL-Multibase / Single	-						
LL-Multibase / Bridge	-						
CMFit Abutment / Ø10	-						
CMFit Abutment / Ø14	CMFit14-NEP-T						
3D-Analog	3D-NEP						
Screw(5ea)	AC-NEP5						

Champions (R)evolution

*set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase			
-		-	-	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Nm) 1.2 Hex			
-				30Nm			
Description-Library code				3501- CHR			
Product	Connection Platform Ø	Ø 3.5/4.5/5.5					
GMS-Scanbody(long)	GMS-CHRL						
LL-Multibase / Single	-						
LL-Multibase / Bridge	-						
CMFit Abutment / Ø10	-						
CMFit Abutment / Ø14	CMFit14-CHR-H						
3D-Analog	3D-CHR						
Screw(5ea)	AC-CHR5						

NEOSS ProActive

Champions (R)evolution

Compatible to Implant Systems

MultiUnit(Multibase)

Dentium Screw Abutment

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
		-		 1.27Hex 10Nm

Description-Library code	0203- MUDE45	0204- MUDE55
Connection Platform Ø	Ø 4.5	Ø 5.5
Product		
SC-Scanbody	SC-MUDE	SC-MUDE
LB-Multibase / Bridge	LB-MUDE45	LB-MUDE55
3D Analog	-	-
Screw(Sea)	AC-MUDE5	AC-MUDE5

Straumann Tissue Level SynOcta coping

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
		-		 Torx 35Nm

Description-Library code	0706- MUIT
Connection Platform Ø	Ø 4.8
Product	
SC-Scanbody	SC-MUIT
LB-Multibase / Bridge	LB-MUIT
3D Analog	-
Screw(Sea)	AC-MUITS

Straumann Bonelevel SRA

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
				 Torx 15Nm

Description-Library code	0707- MUSR35	0708- MUSR46
Connection Platform Ø	Ø 3.5	Ø 4.6
Product		
SC-Scanbody	SC-MUSR35	SC-MUSR46
LB-Multibase / Bridge	LB-MUSR35	LB-MUSR46
3D Analog	3D-MUSR35	3D-MUSR46
Screw(Sea)	AC-MUSR5	AC-MUSR5

Compatible to Implant Systems

MultiUnit (Multibase)

Nobel MultiUnit 4.8

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
				unigrip 15Nm

Description-Library code	0812- MUTB
Connection Platform Ø	Ø4.8
Product	
SC-Scanbody	SC-MUTB / SC2-MUTB
LB-Multibase / Bridge(4.0mm)	LB-MUTB2
LB-Multibase / Bridge(5.5mm)	LB-MUTBL2
3D Analog	3D-MUTB
Screw(5ea)	AC-MUTB5

Ankylos Multi-Unit

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Link	3D-Analog	Screw	Driver (Nm)
		-		1.0Hex 10Nm

Description-Library code	0907- AKLM	0908- AKLR
Connection Platform Ø	Mini	Regular
Product		
SC-Scanbody	SCB-16	SCB-16
LB-Multibase / Bridge	LB-AKLM	LB-AKLR
3D Analog	-	-
Screw(5ea)	AC-AKL5	AC-AKL5

Zimmer Tapered Screw-vent Multi-Unit(TSV)

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Nm)
				1.27Hex 20Nm

Description-Library code	1206- MUTSV
Connection Platform Ø	Universel
Product	
SC-Scanbody	SC-MUTSV
LB-Multibase / Bridge	LB-MUTSV
3D Analog	3D-MUTSV
Screw(5ea)	AC-MUTSV5


MultiUnit (Multibase)

GTS Transfer Scanbody

GTS Transfer Scanbody

Digital implant technology

Convenient & prompt digital implant impression
Ideal subgingiva designing referenced from digital library



Differentiate the Gingival height by appearance



- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co guarantees durability which can be autoclavable and reusable.
- Provides various size with **different gingiva heights**.
– 5mm(GH:2mm), 7mm(GH:4mm), 9mm(GH:6mm)
- Osstem TS and Dentium Super line connection available.
- **Highly precise Geo library** is provided for CAD s/w (3shape & exocad)
- **Special tool** allows to find right position and prevent detachment of screw
- **Customized Abutment Design** with CNC or Premilled abutment is available
- **Zirconia prothesis design** with Geo Link.

Guide Tool



Guide Tool protects the scan body from any damage and helps to find dental implant connection easily.




GTS Scanbody is multifunctional that can be also used as a healing abutment and impression coping.

The groove on the scan body allows users to define different gingiva types and provide a bite index when taking an impression.

GTS Transfer scanbody

Direct Impression → Digital process



GTS Transfer Scanbody List



Guide Tool

1,2 Hex driver(SUM/SURO)
1,27 Hex driver(DER)

DER – Series (Dentium Super-line / Megagen AnyOne / Dentis S-clean / Neo IS)

Diameter GH / Height	4,0mm	4,5mm	5,5mm	6,5mm
2,0mm / 5,0mm	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
4,0mm / 7,0mm	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
6,0mm / 9,0mm	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560
Guide Tool	GTST40	GTST45	GTST55	GTST65

* Laser mark [DER55: Diameter 5,5mm]
Define height and diameter by a different groove on the scan body.

SUM – Series (Osstem Implant TS Mini)

Diameter GH / Height	4,0mm	4,5mm
2,0mm / 5,0mm	GTS-SUM4020	GTS-SUM4520
4,0mm / 7,0mm	GTS-SUM4040	GTS-SUM4540
6,0mm / 9,0mm	GTS-SUM4060	GTS-SUM4560
Guide Tool	GTST40	GTST45

* SUM45: Diameter 4,5mm
Define height and diameter by a different groove on the scan body.

SURO – Series (Osstem Implant TS Standard / Dio UF Regular)

Diameter GH / Height	4,5mm	5,0mm	6,0mm	7,0mm
2,0mm / 5,0mm	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
4,0mm / 7,0mm	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
6,0mm / 9,0mm	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
Guide Tool	GTST45	GTST50	GTST60	GTST70


* SURO50: Diameter 5,0mm
Define height and diameter by a different groove on the scan body.

GTS Transfer scanbody

Sirona Ti-base

Sirona Multibase

Sirona Ti-base



Perfect fit for inCoris ZI meso

Improved design for making emergence profile

Zirconia abutment block for Sirona Multibase



Scan Cap for Sirona Multibase

Description

	SHADE S2	SHADE S3
	S-MEZ0-S2	S-MEZ0-S3

Scan Cap

Scan Cap (short)	Scan Cap (long)
SCSN-CAP	SCSNL-CAP

Ossstem Implant TS

		Platform Ø Connection	TS Mini		Ø3.0/3.5	TS Standard		Ø4.0/4.5/5.0	Driver (Ncm)	
	Single [GH:1.2mm]		SRN2-SUM-H2		SRN2-SURO-H2		20Ncm		Ø3.0/3.5	
	Single [GH:2.6mm]		-		SRN2-SUR026-H2		30Ncm		Ø4.0/4.5/5.0	

*set (Sirona Multibase + screw)

Ossstem Implant US / Biomet 3i Osseotite 4.1

		Platform Ø Connection	US Regular		Ø4.1	Driver (Ncm)	
	Single		SRN2-EXR-H2			30Ncm	

*set (Sirona Multibase + screw)

Sirona Ti-base

Sirona Multibase

Description

Dentium Implantium / SuperLine

Platform Ø Connection	universal	Ø3.7/4.0/4.5	*set (Sirona Multibase + screw)
Single [GH:1.2mm]	SRN2-DER-H2		
Single [GH:2.6mm]	SRN2-DER26-H2		25~30Ncm

Megagen AnyRidge

Platform Ø Connection	Universal	Ø3.5/4.0/4.5	*set (Sirona Multibase + screw)
Single	SRN2-MARR-H2		35Ncm

Megagen AnyOne

Platform Ø Connection	Universal	Ø3.5/4.0/4.5/5.0	*set (Sirona Multibase + screw)
Single [GH:1.2mm]	SRN2-DERA-H2		
Single [GH:2.6mm]	SRN2-DERA26-H2		25Ncm

Neo Biotech IS

Platform Ø Connection		Ø3.5/4.0/4.5	*set (Sirona Multibase + screw)
Single [GH:1.2mm]	SRN2-DERN-H2		
Single [GH:2.6mm]	SRN2-DERN26-H2		30Ncm

Straumann Bone Level

Platform Ø Connection	NC	Ø3.3	RC	Ø4.1/4.8	*set (Sirona Multibase + screw)
Single	SRN2-BLM-R2		SRN2-BLR-R2		35Ncm

Sirona Ti-base

Sirona Multibase

Description

Sirona Ti-base

Straumann Tissue Level

	Platform Ø Connection	RN Ø4.8	WN Ø6.5	*set (Sirona Multibase + screw) Driver (Nm) Torx
Single		SRN2-ITR-02	SRN2-ITW-02	35Nm

Nobel Biocare Active

	Platform Ø Connection	NP Ø3.5	RP Ø4.3/5.0	*set (Sirona Multibase + screw) Driver (Nm) unigrip
Single		SRN2-ACN-H2	SRN2-ACR-H2	35Nm

Nobel Biocare Branemark

	Platform Ø Connection	RP Ø4.1	*set (Sirona Multibase + screw) Driver (Nm) unigrip
Single		SRN2-EXRB-H2	35Nm

Nobel Biocare Replace

	Platform Ø Connection	NP Ø3.5	RP Ø4.3	WP Ø5.0	*set (Sirona Multibase + screw) Driver (Nm) unigrip
Single		SRN2-RPM-E2	SRN2-RPR-E2	SRN2-RPW-E2	35Nm

Dentsply Xive / Frialit2

	Platform Ø Connection	Ø3.4	Ø3.8	Ø4.5	*set (Sirona Multibase + screw) Driver (Nm) 1.22Hex
Single		SRN2-XVM-H2	SRN2-XVR-H2	SRN2-XVW-H2	24Nm

Sirona Ti-base

Sirona Multibase

Description

Dentsply Ankylos

		Platform Ø Connection	universal  Ø 3.4/4.5/5.5	*set (Sirona Multibase + screw) Driver (Ncm)  1.0Hex
Single	Bridge	Single [GH:1.0mm]	SRN2-AKLL1-T	15Ncm
		Single [GH:2.0mm]	SRN2-AKLL2-T	
		Bridge [GH:1.0mm]	SRN2-AKLL1-N	
		Bridge [GH:2.0mm]	SRN2-AKLL2-N	

Sirona Ti-base

Zimmer Tapered Screw-vent

		Platform Ø Connection	 Ø3.5	 Ø4.5	 Ø5.7	*set (Sirona Multibase + screw) Driver (Ncm)  1.25Hex
Single		Single	SRN2-TSV3-H2	SRN2-TSV4-H2	SRN2-TSV5-H2	30Ncm

Astra Tech Osseospeed

		Platform Ø Connection	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0	Driver (Ncm)  1.27Hex
Single		Single	SRN2-ASR-H2	SRN2-ASW-H2	20Ncm 
					25Ncm 

CCM Casting Abutment

CCM Casting Abutment

Product Information

CCM Casting Abutment

Asia's first developed product in 2009
 Six years of accumulated clinical data
 Excellent biocompatibility
 Incl. screw
Optimum platform design for emergence profile
Could make all type of implant restorations



Description

Osstem Implant TS

		Platform Ø Connection		TS Mini Ø3.0/3.5		TS Standard Ø4.0/4.5/5.0		Driver (Ncm) 1.2Hex	
	Single	SUM-H2	GH : 1mm	SURO-H2	20Ncm		SURO22-H2		
			GH : 2.2mm	SURO30-H2			SURO-N2		
			GH : 3mm	SURO22-N2			SURO30-N2		
	Bridge	SUM-N2	GH : 1mm	SURO-N2	30Ncm		SURO22-N2		
			GH : 2.2mm	SURO30-N2			SURO-N2		
			GH : 3mm	SURO22-N2			SURO30-N2		

Osstem Implant US

		Platform Ø Connection		US Mini Ø3.5	US Regular Ø4.1	US Wide Ø5.1	Driver (Ncm) 1.2Hex	
	Single	EXM-H2	EXR-H2	EXW-H2	EXR-N2	EXW-N2		
	Bridge							

Osstem Implant SS

		Platform Ø Connection		SS Regular Ø4.8	SS Wide Ø6.0	Driver (Ncm) 1.2Hex	
	Single	ISR-02	ISW-02	ISR-N2	ISW-N2		
	Bridge						

CCM Casting Abutment

Description

DIO UF / UFII

		Platform Ø Connection	Ø 3.8/4.0/4.5/5.0	*set (CCM Abt + screw)
Single	Bridge		SURD-H2	Driver (Nm) 1.2Hex
			SURD-N2	30Ncm

Neo Biotech IS

		Platform Ø Connection	Ø 3.5/4.0/4.5	*set (CCM Abt + screw)
Single	Bridge		DERN-H2	Driver (Nm) 1.2Hex
			DERN-N2	30Ncm

Straumann Bone Level

		Platform Ø Connection	NC Ø3.3	RC Ø4.1/4.8	*set (CCM Abt + screw)
Single	Bridge		BLM-R2	BLR-R2	Driver (Nm) Torx
			BLM-N2	BLR-N2	35Ncm

Straumann Tissue Level

		Platform Ø Connection	RN Ø4.8	WN Ø6.5	*set (CCM Abt + screw)
Single	Bridge		ITR-02	ITW-02	Driver (Nm) Torx
			ITR-N2	ITW-N2	35Ncm

Nobel Biocare Active

		Platform Ø Connection	NP Ø3.5	RP Ø4.3/5.0	*set (CCM Abt + screw)
Single	Bridge		ACN-H2	ACR-H2	Driver (Nm) unigrip
			ACN-N2	ACR-N2	35Ncm

Nobel Biocare Branemark

		Platform Ø Connection	NP Ø3.5	RP Ø4.1	WP Ø5.1	*set (CCM Abt + screw)
Single	Bridge		EXMB-H2	EXRB-H2	EXWB-H2	Driver (Nm) unigrip
			EXMB-N2	EXRB-N2	EXWB-N2	35Ncm

CCM Casting Abutment

Description

Nobel Biocare Replace

		Platform Ø Connection	NP		Ø3.5	RP		Ø4.3	WP		Ø5.0	Driver (Ncm)	Driver unigrip	*set (CCM Abt + screw)
		Single		RPM-E2		RRP-E2		RPW-E2		RPW-N2				
		Bridge		RPM-N2		RRP-N2		RPW-N2						35Ncm

Nobel Biocare Multi Unit

		Platform Ø Connection			Ø4.8		Driver (Ncm)	Driver unigrip	*set (CCM Abt + screw)
		Bridge		MUTB-N2					15Ncm

Dentsply Xive / Frialit2

		Platform Ø Connection		Ø3.4		Ø3.8		Ø4.5	Driver (Ncm)	Driver 1.22Hex	*set (CCM Abt + screw)
		Single	XVM-H2		XVR-H2		XVW-H2				
		Bridge	XVM-N2		XVR-N2		XVW-N2				24Ncm

Dentsply Ankylos

		Platform Ø Connection		A/B/C/D		Ø3.4/4.5/5.5		Driver (Ncm)	Driver 1.0Hex	*set (CCM Abt + screw)
		Single		-						
		Bridge		DNF3-N						25Ncm

Biomet Bi Certain

		Platform Ø Connection		Ø3.4		Ø4.1		Driver (Ncm)	Driver 1.2Hex	*set (CCM Abt + screw)
		Single	3ICM-H2		3ICR-H2					
		Bridge	3ICM-N2		3ICR-N2					20Ncm

Biomet Bi Osseotite

		Platform Ø Connection		Ø4.1		Ø5.0		Driver (Ncm)	Driver 1.2Hex	*set (CCM Abt + screw)
		Single	EXR-H2		EXT-H2					
		Bridge	EXR-N2		EXT-N2					20Ncm

Zimmer Tapered Screw-vent

		Platform Ø Connection		Ø3.5		Ø4.5		Ø5.7	Driver (Ncm)	Driver 1.25Hex	*set (CCM Abt + screw)
		Single	TSV3-H2		TSV4-H2		TSV5-H2				
		Bridge	TSV3-N2		TSV4-N2		TSV5-N2				30Ncm

CCM Casting Abutment

Description

Zimmer Spline

		Platform Ø Connection	 Ø 3.7 / 4.0	Driver (Nm)
		Single	SP40-H	1.25Hex 30Nm
				

*set (CCM Abt + screw)

Astra Tech Osseospeed

		Platform Ø Connection	X-Small  Ø 3.0	Aqua  Ø 3.5/4.0	Lilac  Ø 4.5/5.0	Driver (Nm)
		Single	AS3-T2	ASR-H2	ASW-H2	1.27Hex 15Nm Ø 3.0
		Bridge	-	ASR-N2	ASW-N2	20Nm Ø 3.5/4.0 25Nm Ø 4.5/5.0
						

*set (CCM Abt + screw)

MIS C1

		Platform Ø Connection	C1  Ø 3.3	C1  Ø 3.75/4.2	C1  Ø 5.0	Driver (Nm)
		Single	C1M-T	C1R-T	C1W-T	1.25Hex 30Nm
		Bridge	C1M-N	C1R-N	C1W-N	
						

*set (CCM Abt + screw)

SIC ace / max

		Platform Ø Connection	 Ø 3.3	 Ø 4.2	Driver (Nm)
		Single	SICM-H	SICR-H	1.2Hex 20Nm
					

*set (CCM Abt + screw)

BEGO

		Platform Ø Connection	 Ø 3.7	 Ø 4.1	Driver (Nm)
		Single	BEG037-H	BEG041-H	1.25Hex 30Nm
		Bridge	BEG037-N	BEG041-N	
					

*set (CCM Abt + screw)

Anthogyr AXIOM

		Platform Ø Connection	universal  Ø 3.4/4.0/4.6/5.2	Driver (Nm)
		Single	ANG-E	1.18Hex 25Nm
		Bridge	ANG-N	
				

*set (CCM Abt + screw)

Tekka In-Kone

		Platform Ø Connection	Standard  Ø 3.5/4.0/4.5/5.5	Driver (Nm)
		Single	TEK-H	1.2Hex 15Nm
				

*set (CCM Abt + screw)

Geo Screw Driver



Angled Driver

Square Type / Ratchet Type

Use with LS Multibase Abutment

Maximum Angle 25°

Recommended torque allowance : below 30Ncm.(reference manufacture's guide)

Caution: unstable connection may bend or destroy the driver

Recommended usage (1.2hex/1.27hex)-10 times

(Universal Type)- 20 times

Connect with Square type Torque wrench

May use Round Type Torque wrench
by connecting HDEC - Short (Clinic)Handle



Connect Straight & Angle driver
with the Torque Wrench



Straight Driver

Square Type / Ratchet Type

Connect by manufacturer's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench
by connecting HDEC - Short (Clinic)Handle

Universal - Round hexagon

Square Type / Ratchet Type

Connect by manufacturer's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench
by connecting HDEC - Short (Clinic)Handle

Angled Driver

3i & Osstem - 1.2mm Hex

ASD3IL Square -20mm	ASD3IR Square -13mm	ALD3IL Ratchet -20mm	ALD3IR Ratchet -13mm

Astra & Zimmer & Dentium - 1.27mm Hex

ASDZIL Square -20mm	ASDZIR Square -13mm	ALDZIL Ratchet -20mm	ALDZIR Ratchet -13mm

Universal - Round hexagon

ASDUVL Square -20mm	ASDUVR Square -13mm	ALDUVL Ratchet -20mm	ALDUVR Ratchet -13mm

Geo Screw Driver

Straight Driver

3i - 1.2mm Hex

			
SSD3IL Square -20mm	SSD3IR Square -13mm	SLD3IL Ratchet -20mm	SLD3IR Ratchet -13mm
	SSD3IEL Square -35mm		

Zimmer - 1.27mm Hex

			
SSDZIL Square -20mm	SSDZIR Square -13mm	SLDZIL Ratchet -20mm	SLDZIR Ratchet -13mm

Astra - 1.27mm Hex

			
SSDATL Square -20mm	SSDATR Square -13mm	SLDATL Ratchet -20mm	SLDATR Ratchet -13mm
	SSDATEL Square -35mm		

Nobel - unigrip

			
SSDNBL Square -20mm	SSDNBR Square -13mm	SLDNBL Ratchet -20mm	SLDNBR Ratchet -13mm
	SSDNBEL Square -35mm		

Straumann - Torx

			
SSDSTL Square -20mm	SSDSTR Square -13mm	SLDSTL Ratchet -20mm	SLDSTR Ratchet -13mm
	SSDSTEL Square -35mm		

Geo Screw Driver

Geo Screw Driver

Driver Tool

HDEL	Long (Lab) Handle		
HDEC	Short (Clinic) Handle		
GTW	Torque Wrench		 Round Type Torque Wrench & Handle Short & Driver  Square Type Torque Wrench & Driver

Lab Driver Set

Description	Product		
	Size	code	Implant System
 1.2 Hex	Blue	Ossstem Implant / Dio Implant Biomet 3i / Xive	
 1.27 Hex	Red	Dentium /Astra Tech Zimmer	
 Unigrip	Yellow	Nobel Biocare	
 Torx	Green	Straumann	
 1.0Hex	Purple	Ankylos	



Lab Driver Set

Instructions of CAD/CAM Abutment

< Procedure for intra-oral scan >



1. Scanbody insert



2. Scanning



3. Matching [Scanbody and Library]



4. Design - 3shape / Exocad



5. Sintering



6. Setting

< Procedure for model scan >





1. Scanbody insert




2. Scanning



3. Matching [Scanbody and Library]




4. Design - 3shape / Exocad




5. Sintering and model setting

Instructions of CCM Casting Abutment

< Screw Type manufacturing >




< Customized Type manufacturing >




Instructions of CCM Casting Abutment




< Hybrid Type manufacturing>



< Bar Overdenture manufacturing>



Geo Library



Free Geo Library available for 3Shape and exocad
High quality of stability and accuracy
One scanbody can be used to design all type of implant restorations
Only implant library including
Ti-base, implant blank, implant bridge



Benefits :

The Geo Library enables you to design individual prosthetic solutions for almost all current implant systems.
The Geo Library is available 3Shape and Exocad libraries, providing 71 different models of titanium bases.
The Geo Library is the dental optimized for CAD was designed as a function.
You are able to fabricate individual milled structures in your own lab, for glued or one piece structures with maximum precision and individuality.



Application :

Individual ceramic abutments on titanium bases

Individual one-piece titanium abutments

Occlusal screwable crowns and bridges

Bars, telescopes

Screw retained abutment

Customized Zirconia abutment

Customized titanium abutment



Selection of Geo Library

Selection of Library

Library for Implant System

- Geo_0101_SUM
- Geo_0102_SURO
- Geo_0103_EXM
- Geo_0104_EXR

Geo product code for different type

Cement gap option with correspond Scanbody

Category	Code	Scanbody Type
Cement [gap]-Zirconia	027	SC
	035	SC
	042	SC
	052	SC
	027	GMS-S
	035	GMS-S
	045	GMS-S
	052	GMS-S
	027	GMS-L
	035	GMS-L
Cement [gap]-PMMA,Metal	042	GMS-L
	052	GMS-L
	08	PMMA-SC
	10	PMMA-SC
	08	PMMA-GMS-S
	10	PMMA-GMS-S
	08	PMMA-GMS-L
	10	PMMA-GMS-L
	08	PMMA-GMS-S
	10	PMMA-GMS-L

Wizard Detect Implant Position

ONE-CLICK REGISTRATION
SELECT IMPLANT FOR TOOTH 21
Geo_0101 SUM
CMFit
10mm SC

LS(short)
LL_C13_H40
LL_C13_H55
LL_C20_H43
LL_C20_H60
CMFit


Click marked point on scan
Best fit matching

Back Next

Back Next

Back Next

Back Next



Geo BS Library- Bridge&Single

GEOMEDI

FOLLOW-ME!
TECHNOLOGY GROUP
hyperDENT
COMPACT · CLASSIC · OPTIONS



Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
Osstem	TS M	SUM	√	√	SC-SUM	SUMC	√
	TS Standard	SURO	√	√	SC-SURO	SUROC	√
	US M	EXM	√	√	SC-EXM	EXMC	√
	US R	EXR	√	√	SC-EXR	EXRC	√
	US W	EXW	√	√	SC-EXW	EXWC	√
Dentium	Superline	DER	√	√	SC-DER	DERC	√
	Screw-Abutment	MUDE45	√		SC-MUDE	MUDEC	√
		MUDE55	√		SC-MUDE	MUDEC	√
Straumann	Bonelevel 3.3	BLM	√	√	SC-BLM	BLMC	√
	Bonelevel 4.1	BLR	√	√	SC-BLR	BLMC	√
	Bonelevel 4.8	BLR	√	√	SC-BLR	BLRC	√
	Tissue level 4.8	ITR	√	√	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	√
	Tissue level 6.5	ITW	√	√	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	√
	SynOcta coping	MUIT	√		SC-MUIT	MUITC	√
	SRA Abutment	MUSR35	√		SC-MUSR35	MUSRC	√
		MUSR46	√		SC-MUSR46	MUSRC	√
Nobel	Multi Unit	MUTB	√		SC-MUTB	MUTBC	√
	Active 3.5	ACN	√	√	SC-ACN	ACNC	√
	Active 4.3	ACR	√	√	SC-ACR	ACRC	√
	Active 5.0	ACR	√	√	SC-ACR	ACRC	√
	Active 5.5	ACW	√	√	SC-ACW	ACRC	√
	Branemark 3.4	EXM	√	√	SC-EXM(B)	EXMBC	√
	Branemark 4.1	EXR	√	√	SC-EXR(B)	EXRBC	√
	Branemark 5.1	EXW	√	√	SC-EXW(B)	EXWBC	√
	Replace 3.5 rpm	RPM	√	√	SC-RPM	RPMC	√
	Replace 4.3 rpm	RPR	√	√	SC-RPR	RPRC	√
	Replace 5.0 rpm	RPW	√	√	SC-RPW	RPRC	√
	Replace 6.0 rpm	RP60	√	√	SC-RP60	RPRC	√
Dentsply	Xive 3.0	XV3	√	√	SC-XV3	XV3C	√
	Xive 3.4	XVM	√	√	SC-XVM	XVC	√
	Xive 3.8	XVR	√	√	SC-XVR	XVC	√
	Xive 4.5	XVW	√	√	SC-XVW	XVC	√
	Xive 5.5	XV55	√	√	SC-XVW	XVC	√
Ankylos	AKII	AKLL	√		OS-AKLL	DNF3C	special screw

Geo BS Library- Bridge&Single

Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	GEO screw code	Compatible with original
3i	Certain 3.4	BICM	✓	✓	SC-BICM	3ICC	v
	Certain 4.1	BICR	✓	✓	SC-BICR	3ICC	v
	Certain 5.0	BICW	✓	✓	SC-BICR	3ICC	v
	Certain 6.0	BIC60	✓	✓	SC-BICR	3ICC	v
	External 3.4	EX3M	✓	✓	SC-EX3M	EXRC	v
	External 4.1	EXR	✓	✓	SC-EXR	EXRC	v
	External 5.0	EXT	✓	✓	SC-EXR	EXRC	v
	External 6.0	EXTW	✓	✓	SC-EXR	EXRC	v
Zimmer	TSV 3.5	TSV3	✓	✓	SC-TSV3	TSVC	v
	TSV 4.5	TSV4	✓	✓	SC-TSV4	TSVC	v
	TSV 5.7	TSV4	✓	✓	SC-TSV5	TSVC	v
Astra	Astra EV 3.0	AEV30	✓		SC-AEV30	AEV30C	v
	Astra EV 3.6	AEV36	✓		SC-AEV36	AEV36C	v
	Astra EV 4.2	AEV42	✓		SC-AEV42	AEV42C	v
	Astra EV 4.8	AEV48	✓		SC-AEV48	AEV4854C	v
	Astra EV 5.4	AEV54	✓		SC-AEV54	AEV4854C	v
	Astra 3.5	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.0	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.5	ASW	✓	✓	SC-ASW	ASWC	v
	Astra 5.0	ASW	✓	✓	SC-ASW	ASWC	v
Camlog	Camlog 3.3	CAM33	✓		SC-CAM33	CAMC	v
	Camlog 3.8	CAM38	✓		SC-CAM38	CAMC	v
	Camlog 4.3	CAM43	✓		SC-CAM43	CAMC	v
	Camlog 5.0	CAM50	✓		SC-CAM5060	CAM5060C	v
	Camlog 6.0	CAM60	✓		SC-CAM5060	CAM5060C	v
Magazen	AnyOne	MAO	✓	✓	SC-DER	MAOC	v
	Anyridge	MARR	✓	✓	SC-MARR	MARRC	v
Dentis	S-clean	DET	✓	✓	SC-DER	DERC	v
Dio	UF	SURD	✓	✓	SC-SURO	SURDC	v
Neo	IS	DERN	✓	✓	SC-SURO	DERNC	v
MIS	C1 3.3	C1M	✓		SC-C1M	C1MC	v
	C1 3.75 / 4.2	C1R	✓		SC-C1R	C1C	v
	C1 5.0	C1W	✓		SC-C1W	C1C	v

Geo Product List (Scanbody / Multibase / CMFit / 3D Analog / CCM Abutment / Strona Ti-base)

No	Company	Implant system	Platform	Library Code	Type	sc(os) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Strona Ti-base
101			3.0/3.5	0101-SUM	Single	SC-SUM	GMS-SUMS	LL2-SUM-M13-H2	LS2-SUM-MH2	CMFIT10-SUM-H2	CMFIT14-SUM-H2	3-D-SUM	AC-SUM5	SUM-H2	SRN2-SUM-H2
102	TS	4.0/4.5/5.0-	0102-SURO	Single	Bridge	GMS-SUML	LL2-SUM20-H2	LS2-SUM-N2	CMFIT10-SUM-N2	-	-	-	-	SUM-N2	-
103	Oststem	3.5	0103-EXM	Single	Bridge	GMS-SUROS	LL2-SUM13-N2	LS2-SUM-N2	CMFIT10-SURO-H2	CMFIT14-SURO-H2	3D-SURO	AC-SUROS	SURO-H2	SRN2-SURO-H2	
104	US	4.1	0104-EXR	Single	Bridge	GMS-SUROL	LL2-SURO20-H2	LS2-SURO-N2	CMFIT10-SURO-N2	CMFIT14-SURO-N2	3D-SURO	AC-SUROS	SURO30-H2	SRN2-SURO30-H2	
105		5.1	0105-EXW	Single	Bridge	GMS-EXWS	LL2-SURO13-N2	LS2-SURO-N2	CMFIT10-SURO-N2	CMFIT14-SURO-N2	3D-SURO	AC-SUROS	SURO22-N2	-	
106		4.8	0106-4SR	Single	Bridge	GMS-EXML	LL2-EXW-H2	LS2-EXW-H2	CMFIT10-EXW-H2	CMFIT14-EXW-H2	3D-EXW	AC-EXM5	EXW-H2	SRN2-EXW-H2	
107	SS	6.0	0107-4SW	Single	Bridge	GMS-EXRS	LL2-EXR-H2	LS2-EXR-H2	CMFIT10-EXR-H2	CMFIT14-EXR-H2	3D-EXR	AC-EXR5	EXR-H2	SRN2-EXR-H2	
201	Dentium	3.7/4.0/4.5 (Universal)	0201-DER	Single	Bridge	GMS-DEWS	LL2-EXW-N2	LS2-EXW-N2	CMFIT10-EXW-N2	CMFIT14-EXW-N2	3D-EXW	AC-EXW5	EXW-N2	-	
202	NR-Line	3.1/3.6/4.3 (Standard)	0202-DNR	Single	Os-DNR	GMS-MARRS	LL2-EXR-N2	LS2-EXR-N2	CMFIT10-EXR-N2	CMFIT14-EXR-N2	3D-DER	AC-DERS	DER-H2	SRN2-DER-H2	
301	AnyRidge	3.5/4.0/4.5 (Universal)	0301-MARR	Single	Bridge	GMS-MARR	LL2-MARR20-H2	LS2-MARR-H2	CMFIT10-MARR-H2	CMFIT14-MARR-H2	3D-MARR	AC-MARR5	MARR-H2	SRN2-4MARR-H2	
302	Megagen	3.0/4.0/4.5/5.0 (Universal)	0302-DERA	Single	SO-DER	GMS-DERS	LL2-MARR20-N2	LS2-MARR-N2	-	-	-	-	MARR-N2	-	
303	EZ-Plus	3.3	0303-MEPM	Single	Bridge	GMS-DERL	LL2-DERA13-H2	LS2-DERA-N2	CMFIT10-DERA-H2	CMFIT14-DERA-H2	3D-DER	AC-MAR5	DER-H2	SRN2-DERA-H2	
304	EX-feel	4.0/4.5/5.0 (Regular)	0304-MEPR	Single	SC-MEPM	GMS-MEPR	LL2-DERA20-H2	LS2-DERA20-N2	CMFIT10-MEPM-E	CMFIT14-MEPM-E	-	AC-MEPMS	MEPM-E	-	
305		4.8	0305-ITRK	Single	Bridge	SC-ITRW	LL2-ITRK<20	LS2-ITRK<20	CMFIT10-ITRK<20	CMFIT14-ITRK<20	3D-ITR	AC-EXR5	ITRK-C02	-	
306		Mini	3.0/3.25	0306-MMN	Single	SC-MMN	-	-	CMFIT10-MMN-H	-	-	AC-MMNS	-	-	
401	Dio	SM	3.8/4.1	0401-DSMM	Single	SC-DSMM	-	-	CMFIT10-DSMM-T	-	-	AC-DSMM5	DSMM-T	-	
402		4.5/5.0	0402-DSMR	Single	SC-DSMR	GMS-DSMR	-	-	CMFIT10-DSMR-T	-	-	AC-DSMR5	DSMM-N	-	

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
403	Dio	3.0/3.3 (Narrow)	0403-DN30	Single	SC-MMN	-	-	-	-	CMFIT10-DN30-H	-	-	-	AC-DN305	-
404	UF / UFII	3.8/4.0/4.5/5.0 (Standard)	0404-SURD	Single	SC-SURO	GMS-SUROS	LL2-SURD-H2	CMFIT10-SURD-H2	CMFIT14-SURD-H2	-	-	SURD+H2	-	-	-
501	NeoBiotech	IS	0501-DERN	Single	SC-SURO	GMS-DERNS	LL2-SURD20-N2	-	-	CMFIT10-SURD-N2	CMFIT14-SURD-N2	3D-SURD	AC-SURD5	SURD+N2	-
601	Dentis	3.5/4.0/4.5/5.0-	0601-DERT	Single	SC-DER	GMS-DERS	LL2-DERNH2	-	-	CMFIT10-DERNH2	CMFIT14-DERNH2	3D-DER	AC-DERNS	DERNH2	SRN2-DERNH2
709		2.9 (SC)	0700-BLS	Single	-	GMS-BLSL	LL2-BLS-R	-	-	CMFIT10-BLS-R	CMFIT14-BLS-R	3D-BLS	AC-BLS5	-	-
701	Straumann	BoneLevel	3.3 (NC)	0701-BLM	Single	SC-BLM	GMS-BLMS	LL2-BLM-N2	LS2-BLM-R2	CMFIT10-BLM-R2	CMFIT14-BLM-R2	3D-BLM	AC-BLM5	BL-M-R2	SRN2-BLM-R2
702		4.1/4.8 (RC)	0702-BLR	Single	SC-BLR	GMS-BLRL	LL2-BLR-N2	LL2-BLR20-R2	LS2-BLR-R2	CMFIT10-BLR-R2	CMFIT14-BLR-R2	3D-BLR	AC-BLR5	BLR-R2	SRN2-BLR-R2
703		3.5 (NN)	0703-IT35	Single	SC-IT35	-	-	-	LS2-IT35-O	CMFIT10-IT35-O	-	-	AC-IT35	-	-
704	Straumann	Tissue Level	4.8 (RN)	0704-ITR	Single	GMS-ITRS	LL2-ITR20-O2	LS2-ITR-O2	CMFIT10-ITR-O2	CMFIT14-ITR-O2	3D-ITR	-	-	ITR-O2	SRN2-ITR-O2
705		6.5 (WN)	0705-ITW	Single	SC-ITRW	GMS-ITRW	LL2-ITR20-N2	LL2-ITR-N2	-	-	-	-	ITR-N2	-	-
710		BLX	3.6 - 6.5	0710-BLX	Single	GMS-BLXL	LL2-BLX-T	-	-	CMFIT10-BLX-T	CMFIT14-BLX-T	3D-BLX	AC-BLX5	-	-
801	Nobel Biocare	3.0	0801-AC30	Single	OS-AC30	-	LL2-AC30-T2	LS2-AC30-T2	CMFIT10-AC30-T2	CMFIT14-AC30-T2	3D-AC30	AC-AC305	-	-	
802		3.5 (NP)	0802-ACN	Single	SC-ACN	GMS-ACNS	LL2-ACNH2	LS2-ACNH2	CMFIT10-ACNH2	CMFIT14-ACNH2	3D-ACN	AC-ACN5	ACN-H2	SRN2-ACN-H2	
803		4.3/5.0 (RP)	0803-ACR	Bridge	SC-ACR	GMS-ACRS	LL2-ACRN2	LS2-ACRN2	CMFIT10-ACR-H2	CMFIT14-ACR-H2	3D-ACR	AC-ACR5	ACR-H2	SRN2-ACR-H2	
804		5.5 (WP)	0804-ACW	Single	SC-ACW	GMS-ACW	LL2-ACNW2	LS2-ACNW2	CMFIT10-ACNW2	CMFIT14-ACNW2	3D-ACW	-	-	-	
805	Branemark	3.5 (NP)	0805-EXMB	Single	SC-EXMB	-	LL2-EXMB-H2	LS2-EXMB-H2	CMFIT10-EXMB-H2	CMFIT14-EXMB-H2	3D-EXM	AC-EXMB5	EXMB-H2	-	
806		4.1 (RP)	0806-EXRB	Bridge	SC-EXRB	-	LL2-EXRB-H2	LS2-EXRB-H2	CMFIT10-EXRB-H2	CMFIT14-EXRB-H2	3D-EXR	AC-EXRB5	EXRB-H2	SRN2-EXRB-H2	
807		5.1 (WP)	0807-EXWB	Bridge	SC-EXWB	-	LL2-EXWB-H2	LS2-EXWB-H2	CMFIT10-EXWB-H2	CMFIT14-EXWB-H2	3D-EXW	AC-EXWB5	EXWB-H2	-	

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFIT 10mm	CMFIT 14mm	3D Analog	Screw(Sea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
808			3.5 (NP)	0808-RPM	Single	GMS-RPMS	LL2-RPM-E2	CMFIT10-RPM-E2	CMFIT14-RPM-E2	3D-RPM	AC-RPM5	RPM-E2		SRN2-RPM-E2	
809	Nobel Biocare	Replace Select	4.3 (RP)	0809-RPR	Bridge	SC-RPM	GMS-RPML	LL2-RPM-N2	LS2-RPM-E2	CMFIT10-RPR-E2	CMFIT14-RPR-E2	3D-RPR		RPR-E2	SRN2-RPR-E2
810			5 (WP)	0810-RPW	Single	SC-RPR	GMS-RPPRS	LL2-RPR-E2	LS2-RPR-E2	CMFIT10-RPR-E2	CMFIT14-RPR-E2	3D-RPR		RPR-N2	
811			6.0	0811-RP60	Bridge	SC-RPW	GMS-RPPWS	LL2-RPW-E2	LS2-RPW-E2	CMFIT10-RPW-E2	CMFIT14-RPW-E2	3D-RPW		RPW-E2	SRN2-RPW-E2
901			3.0	0901-XV3	Single	SC-RP60	GMS-RP60S	LL2-RP60-E2	LS2-RP60-E2	CMFIT10-RP60-E2	CMFIT14-RP60-E2	3D-RP60		RPW-N2	
902	Xive		3.4	0902-XVM	Bridge	SC-XVM	-	LL2-XVM-N2	LS2-XVM-H2	CMFIT10-XVM-H2	CMFIT14-XVM-H2	3D-XVM		XVM-H2	SRN2-XVM-H2
903			3.8	0903-XVR	Bridge	SC-XVR	-	LL2-XVR-N2	LS2-XVR-H2	CMFIT10-XVR-H2	CMFIT14-XVR-H2	3D-XVR		XVR-H2	SRN2-XVR-H2
904	Dentsply		4.5	0904-XW	Bridge	SC-XW	-	LL2-XW-N2	LS2-XW-H2	CMFIT10-XW-H2	CMFIT14-XW-H2	3D-XW		XW-H2	SRN2-XW-H2
905			5.5	0905-XW55	Bridge	SC-XW	-	LL2-XW55-N2	LS2-XW55-H2	CMFIT10-XW55-H2	CMFIT14-XW55-H2	3D-XW5		XW55-H2	SRN2-XW55-H2
906	Antylos	3.4/4.5/5.5 (Universal)	0906-AKLL		Single	OS-AKLL	-	LL2-AKLL-T	CMFIT10-AKLL-T	-	-			SRN2-AKLL1-T	
1001	Certain		3.4	1001-3ICM	Bridge	SC-3CM	GMS-3ICMS	LL2-3ICM-N2	LS2-3ICM-H2	CMFIT10-3ICM-H2	-	3D-3ICM		3ICM-H2	
1002			4.1	1002-3ICR	Bridge	SC-3CR	GMS-3ICRL	LL2-3ICR-N2	LS2-3ICR-H2	CMFIT10-3ICR-H2	CMFIT14-3ICR-H2	3D-3ICR		3ICR-H2	
1003			5.0	1003-3ICW	Bridge	SC-3ICR	GMS-3ICRS	LL2-3ICW-N2	LS2-3ICW-H2	CMFIT10-3ICW-H2	CMFIT14-3ICW-H2	3D-3ICW		3ICR-N2	
1004	Biomet 3i		6.0	1004-3IC60	Bridge	SC-3ICR	GMS-3ICRL	LL2-3IC60-N2	LS2-3IC60-H2	CMFIT10-3IC60-H2	CMFIT14-3IC60-H2	3D-3IC60			
1005			3.4	1005-EX3M	Bridge	SC-EX3M	-	LL2-EX3M-H2	LS2-EX3M-H2	CMFIT10-EX3M-H2	-	-	-		
1006	External / Osseotite		4.1	1006-EXR	Bridge	SC-EXR	GMS-EXRS	LL2-EXR-H2	LS2-EXR-H2	CMFIT10-EXR-H2	CMFIT14-EXR-H2	3D-EXR		EXR-H2	SRN2-EXR-H2
1007			5.0	1007-EXT	Bridge	SC-EXR	GMS-EXRS	LL2-EXT-H2	LS2-EXT-H2	CMFIT10-EXT-H2	CMFIT14-EXT-H2	-		EXT-H2	SRN2-EXT-H2
1008			6.0	1008-EXTW	Bridge	SC-EXTW	GMS-EXRS	LL2-EXTW-H2	-	CMFIT10-EXTW-H2	CMFIT14-EXTW-H2	-		EXT-N2	
1101	Implant Direct	Legacy	3.0	1101-TSM	Single	OS-TSM	-	LL2-TSM-H2	-	CMFIT10-TSM-H2	-	3D-TSM	AC-TSM5	-	
1201	Zimmer	Tapered Screw-Vent	3.5	1201-TS/3	Bridge	SC-TS/3	GMS-TS/3	LL2-TS/3-H2	LS2-TS/3-H2	CMFIT10-TS/3-H2	CMFIT14-TS/3-H2	3D-TS/3		TS/3-H2	SRN2-TS/3-H2
1202			4.5	1202-TS/4	Bridge	SC-TS/4	GMS-TS/4L	LL2-TS/4-H2	LS2-TS/4-H2	CMFIT10-TS/4-H2	CMFIT14-TS/4-H2	3D-TS/4		TS/4-N2	SRN2-TS/4-H2
1203-B 11102; Q150			5.7	1203-TS/5	Bridge	SC-TS/5	GMS-TS/5S	LL2-TS/5-H2	LS2-TS/5-H2	CMFIT10-TS/5-H2	CMFIT14-TS/5-H2	3D-TS/5		TS/5-H2	SRN2-TS/5-H2

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	L1-Link	Ls-Link	Cmf1t 10mm	3D Analog	Screw(Sea)	CCM Abutment (Co-Cr base)	Sirona Ti-base
1204	Zimmer	Spline	3.7/4.0	1204-SP40	Single	SC-SP40	-	L1-2-SP40-H	-	-	Cmf1t 4-SP40-H	-	SP40-H	-
1205			5.0	1205-SP50	Single	-	-	L1-2-SP50-H	-	-	Cmf1t 4-SP50-H	-	A-C-TSv6	-
1301			3.0 (Yellow)	1301-A33	Single	SC-A33	-	L1-2-AS3-T2	-	-	Cmf1t 10-AS3-T2	-	A-C-AS35	A33-T2
1302	Osseospeed		3.5/4.0 (Aqua)	1302-ASR	Single	SC-ASR	GMS-ASRS	L1-2-ASR2-H2	L2-ASR2-H2	Cmf1t 10-ASR-H2	Cmf1t 4-ASR-H2	3D-ASR	A-C-ASR5	ASR-H2
1303	Astra Tech		4.5/5.0 (Lilac)	1303-ASW	Bridge	GMS-ASRL	GMS-ASWS	L1-2-ASR2-N2	L2-ASR2-N2	Cmf1t 10-ASR-H2	Cmf1t 4-ASW-H2	3D-ASW	A-C-ASW5	SRN2-ASR-H2
1304			3 (Green)	1304-AE/30	Single	SC-ASW	GMS-ASWL	L1-2-ASW-H2	L2-ASW-H2	Cmf1t 10-ASW-H2	Cmf1t 4-ASW-H2	3D-ASW	A-C-ASW5	ASW-N2
1305			3.6 (Purple)	1305-AE/36	Bridge	SC-AEV30	-	L1-2-AEV30-N2	-	Cmf1t 10-AE/30-T2	-	3D-AEV30	AC-AEV305	-
1306			4.2 (Yellow)	1306-AEV42	Single	SC-AEV42	-	L1-2-AEV30-T2	-	Cmf1t 10-AE/36-T2	Cmf1t 4-AEV36-T2	3D-AEV36	AC-AEV365	-
1307			4.8 (Blue)	1307-AEV48	Bridge	SC-AEV48	-	L1-2-AEV42-T2	-	Cmf1t 10-AE/42-T2	Cmf1t 4-AEV42-T2	3D-AEV42	AC-AEV425	-
1308			5.4 (Brown)	1308-AEV54	Bridge	SC-AEV54	-	L1-2-AEV42-N2	-	Cmf1t 10-AE/42-N2	Cmf1t 4-AEV48-T2	3D-AEV48	AC-AEV48545	-
1401			3.3	1401-CAM33	Single	SC-CAM33	-	L1-2-CAM33-T2	-	Cmf1t 10-AE/48-T2	Cmf1t 4-AEV54-T2	3D-AEV54	-	-
1402			3.8	1402-CAM38	Single	SC-CAM38	-	L1-2-CAM38-T2	-	Cmf1t 10-CAM38-T2	Cmf1t 4-CAM38-T2	3D-CAM38	AC-CAM5	-
1403	Camlog		4.3	1403-CAM43	Bridge	SC-CAM43	-	L1-2-CAM43-T2	-	Cmf1t 10-CAM43-T2	Cmf1t 4-CAM43-T2	3D-CAM43	-	-
1404	Camlog		5.0	1404-CAM50	Single	SC-CAM50	-	L1-2-CAM50-T2	-	Cmf1t 10-CAM50-T2	Cmf1t 4-CAM50-T2	3D-CAM50	-	-
1405			6.0	1405-CAM60	Bridge	SC-CAM560	-	L1-2-CAM60-T2	-	Cmf1t 10-CAM60-T2	Cmf1t 4-CAM60-T2	3D-CAM560	AC-CAMS5605	-
1406			3.3	1406-CN33	Single	GMS-CN33L	-	-	Cmf1t 10-CN33-E	-	3D-CN33	-	-	-
1407			3.8/4.3	1407-CN843	Single	-	GMS-CN843L	-	-	Cmf1t 14-CN843-E	3D-CN843	AC-CN5	-	-
1408			5.0	1408-CN50	Single	-	GMS-CN50L	-	-	Cmf1t 10-CN50-E	3D-CN50	AC-CN505	-	-
1501			3.3 (Narrow)	1501-MS33	Bridge	SC-MS33	-	L1-2-MS33-H	L2-2-MS33-H	Cmf1t 10-MS33-H	-	3D-MS33	AC-MS335	-
1502	Seven	3.7/4.2 (Standard)	1502-MSR	Single	SC-TSV3	GMS-TSV3	L1-2-MSR-H	Cmf1t 10-MSR-H	Cmf1t 14-MSR-H	Cmf1t 14-MSR-H	3D-TSv3	-	-	-
1503	MIS	5.0/5.0 (Wide)	1503-MSW	Single	SC-TSV4	GMS-TSV4	L1-2-MSR-N	Cmf1t 10-MSV-H	Cmf1t 14-MSV-H	Cmf1t 14-MSV-H	3D-TSv4	AC-MSS	-	-
1504		3.3	1504-C1M	Single	SC-C1M	GMS-C1MS	L1-2-C1M-T	Cmf1t 10-C1M-T	Cmf1t 14-C1M-T	Cmf1t 14-C1M-T	3D-C1M	A-C-C1M5	C1M-T	C1M-N
1505		3.75/4.2	1505-C1R	Single	SC-C1R	GMS-C1RS	L1-2-C1R-T	Cmf1t 10-C1R-T	Cmf1t 14-C1R-T	Cmf1t 14-C1R-T	3D-C1R	C1R-T	C1R-N	C1W-T
1506		5.0	1509-C1W	Single	SC-C1W	GMS-C1WS	L1-2-C1W-T	Cmf1t 10-C1W-T	Cmf1t 14-C1W-T	Cmf1t 14-C1W-T	3D-C1W	A-C-C15	C1W-T	C1W-N

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(5ea)	CCM Abutment (Co-Cr base)	Strona Ti-base
1601	SIC	ace / SIC max	3.3	1601-SICM	Single	-	-	LL2-SICM-H	-	CMFIT10-SICM-H	CMFIT14-SICM-H	3D-SICM	SICM-H	SICM-H	-
1602			4.2	1602-SICR	Bridge	SC-SIC	-	LL2-SICR-N	-	CMFIT10-SICR-H	CMFIT14-SICR-H	3D-SICR	SICR-H	SICR-H	-
1701	Thommen	SPI	3.5	1701-SPI35	Bridge	SC-SPI35	-	LL2-SPI35-H	LS2-SPI35-H	CMFIT10-SPI35-H	-	3D-SPI35	AC-SPI35	-	-
1702			4.0	1702-SPI40	Bridge	SC-SPI40	-	LL2-SPI40-H	LS2-SPI40-H	CMFIT10-SPI40-H	CMFIT14-SPI40-H	3D-SPI40	-	-	-
1703	Thommen	SPI	4.5	1703-SPI45	Bridge	SC-SPI45	-	LL2-SPI45-H	LS2-SPI45-H	CMFIT10-SPI45-H	CMFIT14-SPI45-H	3D-SPI45	-	-	-
1704			5.0	1704-SPI50	Bridge	SC-SPI50	-	LL2-SPI50-H	LS2-SPI50-H	CMFIT10-SPI50-H	CMFIT14-SPI50-H	3D-SPI50	AC-SPI5	-	-
1705	Bego	Bego	6.0	1705-SPI60	Bridge	SC-SPI60	-	LL2-SPI60-H	LS2-SPI60-H	CMFIT10-SPI60-H	CMFIT14-SPI60-H	-	-	-	-
1801			3.7	1801-BEGO37	Bridge	SC-BEGO37	-	LL2-BEGO37-H	LS2-BEGO37-H	CMFIT10-BEGO37-H	-	3D-BEGO37	BEGO37-H	-	-
1802	Bego	Bego	4.1	1802-BEGO41	Bridge	SC-BEGO	-	LL2-BEGO41-H	LS2-BEGO41-H	CMFIT10-BEGO41-H	-	-	BEGO41-H	-	-
1803			4.5	1803-BEGO45	Single	SC-BEGO45	-	LL2-BEGO45-H	-	CMFIT10-BEGO45-H	-	-	AC-BEGO5	-	-
1804	Mytis	Arrow	5.5	1804-BEGO55	Single	SC-BEGO55	-	LL2-BEGO55-H	-	CMFIT10-BEGO55-H	-	-	-	-	-
1901			3.74.0	1901-AR3740	Single	OS-AR3740	-	-	-	CMFIT10-AR3740-H	-	-	AC-AR5	-	-
1902	Biohorizon	POI	4.7/6.0	1902-AR4760	Single	OS-AR4760	-	-	-	CMFIT10-AR4760-H	-	-	AC-POI375	-	-
2001			3.7	2001-POI37	Single	SC-POI37	-	LL2-POI37-H	-	CMFIT10-POI37-H	-	-	-	-	-
2002	KYOCERA	Laser Lok	4.2	2002-POI42	Single	SC-POI42	-	LL2-POI42-H	-	CMFIT10-POI42-H	-	-	AC-POI425	-	-
2101			3.0	2101-BH30	Single	OS-BH30	-	LL2-BH30-H	-	CMFIT10-BH30-H	-	-	AC-BH305	-	-
2201	Anthonyar	AXION Bone Level	3.4/4.0/4.6/5.2 (Universal)	2201-ANG	Single	SC-ANG	-	LL2-ANG-E	-	CMFIT10-ANG-E	-	3D-ANG	AC-ANG5	ANG-E	-
2301			4.3/4.8/5.3 (Standard)	2301-WOP	Bridge	SC-WOP	-	-	-	CMFIT10-WOP-O	-	-	AC-WOP5	-	-
2401	Teeka	In-Kone	3.5/4.0/4.5/5.5 (Standard)	2401-TEK	Bridge	SC-TEK	-	-	LS2-TEK-H	CMFIT10-TEK-H	-	3D-TEK	AC-TEK5	TEK-H	-
2501			3.5	2501-KEP35	Bridge	SC-KEP35	-	LL2-KEP35-T	LS2-KEP35-T (SLKEP35)	CMFIT10-KEP35-T	CMFIT14-KEP35-T	3D-KEP35	-	-	-
2502	Keystone	PrimaComex / Max-Tilobe	4.1	2502-KEP41	Bridge	SC-KEP41	-	LL2-KEP41-T	LS2-KEP41-T (SLKEP41)	CMFIT10-KEP41-T	CMFIT14-KEP41-T	3D-KEP41	-	-	-
2503			5.0/6.0	2503-KEP50	Bridge	SC-KEP50	-	LL2-KEP50-T	LS2-KEP50-T (SLKEP50)	CMFIT10-KEP50-T	CMFIT14-KEP50-T	3D-KEP50	AC-KEP5	-	-
2504	Neodent	CM	7.0	2504-KEX50	Bridge	SC-NCM	-	LL2-KEX50-T	LS2-KEX50-T (SLKEX50)	-	CMFIT14-KEX50-T	3D-KEX50	-	-	-
2505			8.0	2505-KEX65	Bridge	SC-NCM	-	LL2-KEX65-T	LS2-KEX65-T (SLKEX65)	-	CMFIT14-KEX65-T	3D-KEX65	-	-	-
2601	Neodent	GM	3.5/4.3/5.0 (Standard)	2601-NCM	Single	SC-NCG	-	LL2-NCG-H	-	CMFIT10-NCM-H	-	3D-NCM	-	-	-
2602			3.5/3.7/4.4/3.5/5.0 (Standard)	2602-NCG	Single	SC-NCG	-	LL2-NCG-H	-	CMFIT10-NCG-H	-	3D-NCG	AC-NCG5	-	-

No	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Link	LS-Link	CMFit 10mm	CMFit 14mm	3D Analog	Screw(See)	CCM Abutment (Co-Cr base)	Sirona TiBase
2701	BICON		2.0	2701-B2C20	Bridge	SC-BC20	-	LL-2-B2C20-N	-	-	-	-	-	-	-
2702	BICON	BICON Implant	2.5	2702-B2C25	Bridge	SC-BC25	-	LL-2-B2C25-N	-	-	-	-	-	-	-
2703			3.0	2703-B2C30	Bridge	SC-BC30	-	LL-2-B2C30-N	-	-	-	-	-	-	-
2801	Medenika	IPS	3.5/4.0/4.5/5.0	2801-MTR	Single	SC-MTR	-	-	-	-	-	-	-	-	AC-MTR5
2901	GC	Adeva Implant (Regular / Wide)	4.0/5.0	2901-GAD	Single	SC-GAD	-	-	-	-	-	-	-	-	AC-GAD5
3001	BIOTECH	KONTACT	3.0	3001-BK30	Single	SC-BK30	-	LL-2-BK30-T	LS2-BK30-T	-	-	3D-BK30	-	-	-
3002			3.6/4.2/4.8/5.4	3002-BK40	Single	SC-BK40	-	LL-2-BK40-T	LS2-BK40-T	-	-	3D-BK40	-	-	-
3101	Medentis medical	ICX	3.75/4.14/4.8	3101-ICX	Single	SC-ICX	-	LL-2-ICX-H	LS2-ICX-H	CMFit14-ICX-H	CMFit14-ICX-H	3D-ICX	AC-ICX5	-	-
3201	B&B Dental	Dura-vit	4.0/3.5/5.0	3201-BNB	Single	-	GMS-BNBL	LL-2-BNB-H	-	-	-	3D-BNB	AC-BNBS	-	-
3301	Bredent Medical	Sky	3.5 / 4.0	3301-SKR	Single	-	GMS-SKRL	LL-2-SKR-T	-	CMFit14-SKR-T	CMFit14-SKR-T	3D-SKR	AC-SKR5	-	-
3401	Neoss	ProActive	3.5/4.0/4.5/5.0/5.5	3401-NEP	Single	-	GMS-NEPL	-	-	-	CMFit14-NEP-T	3D-NEP	AC-NEP5	-	-
3501	Champions	(Revolution	3.5/4.0/4.5/5.5	3501-CHR	Single	-	GMS-CHRL	-	-	CMFit14-CHR-H	CMFit14-CHR-H	3D-CHR	AC-CHR5	-	-

Geo MultiUnit List (Scanbody / Cylinder / 3D Analog)

No	Compatible to implant Systems	Platform	Library Code	Type	SC(OS) Scanbody	LB-short	LB-long	Temporary Cylinder	3D Analog	Screw(See)	CCM Abutment (Co-Cr base)
203	Dentium Screw Abutment	4.5	0203-MUDE5	Bridge	SC-MUDE	LB-MUDE45	LB-MUDE55				AC-MUDES
204		5.5	0204-MUDE5	Bridge	SC-MUIT	LB-MUIT					AC-MUITS
706	Straumann Tissue Level Systech coping	4.8	0706-MUIT	Bridge	SC-MUITS	LB-MUITS35	LB-MUITS45	3D-MUITS35	3D-MUITS46	TA-MUITS46	AC-MUSRS
707	Straumann bonelevel RA	3.4	0707-MUSR35	Bridge	SC-MUSR35	LB-MUSR35	LB-MUSR46	TA-MUSR46	3D-MUSR35	3D-MUSR46	AC-MUSR5
708	Nobel Multi Unit	4.6	0708-MUSR46	Bridge	SC-MUSR46	LB-MUTB2	LB-MUTB2	TA-MUTB	3D-MUTB	AC-MUTBS	MUTB-H
812	Dentsply mini	4.8	0812-MUTB	Bridge	SC-MUTB	SC2-G6N	SC2-G6N				
907	Anteklos Multi-Unit regular	mini	0907-AKLM	Bridge	SC-B16	LB-AKLM	LB-AKLM				AC-AKL5
908	Zimmer Tapered Screw-Vent Tapered Coping	regular	0908-AKLR	Bridge	SC-MUTSV	LB-MUTSV	LB-MUTSV				AC-MUTSV5
1206		4.5	1206-MUTSV	Bridge	SC-MUTSV				3D-MUTSV		

Geo Transfer Scanbody List

No	Compatible to Implant Systems	Platform	Type	Section				
				Diameter GH / Height	4.0mm	4.5mm		
101	Osstem Implant TS Mini	3.0/3.5	Single	2.0 / 5.0	GTS-SUM4020	GTS-SUM4520		
				4.0 / 7.0	GTS-SUM4040	GTS-SUM4540		
				6.0 / 9.0	GTS-SUM4060	GTS-SUM4560		
102	Osstem Implant TS Standard	4.0/4.5/5.0-	Single	Diameter GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
				2.0 / 5.0	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
				4.0 / 7.0	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
201	Dentium Super-Line Implantum	3.7/4.0/4.5-	Single	Diameter GH / Height	4.0mm	4.5mm	5.5mm	6.5mm
				2.0 / 5.0	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
				4.0 / 7.0	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
				6.0 / 9.0	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560

Guide Tool	Diameter	4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.0mm
	code	GTST40	GTST45	GTST50	GTST55	GTST60	GTST65	GTST70

GEOMEDI

Product Catalog - 2022

GeoMedi Co., Ltd.

GEOMEDI

PREMIUM CAD/CAM PRODUCTS

(주) 지 오 메 디

경기도 의왕시 이미로 40 (포일동, 인덕원IT밸리) D동 702~705, 810-2, 811호

Tel : 02-588-2871 / Fax : 02-6280-2849

sales@geomedi.com www.geomedi.com

GeoMedi Co., Ltd.

#811, 810-2, 702-705, Unit D, Indeokwon IT Valley, 40, Imi-ro, Uiwang-si, Gyeonggi-do, Korea

Tel : 82-2-588-2871 Fax : 82-2-6280-2849

sales@geomedi.com www.geomedi.com



GTS Transfer Scanbody

Multi-screw Abutment